

```
In[3]:= << QMRITools`
```

```
In[9]:= Column@QMRIToolsPackages[]
```

```
CardiacTools  
CoilTools  
DenoiseTools  
DixonTools  
ElastixTools  
GeneralTools  
GradientTools  
ImportTools  
IVIMTools
```

```
Out[9]= JcouplingTools  
MaskingTools  
NiftiTools  
PhysiologyTools  
PlottingTools  
ProcessingTools  
RelaxometryTools  
SimulationTools  
TensorTools  
VisteTools
```

```
In[10]:= QMRIToolsFunctions[60]
```

Functions

ADCCalc	DTItoolExpFile	ImportBvalvec	PlotData3D
AddNoise	DTItoolExpInd	ImportBvec	PlotDefGrid
AlignRespLog	DTItoolExpTens	ImportDTI	PlotDuty
AngleCalc	ECalc	ImportExploreDTITens	PlotIVIM
AngleMap	ECVCalc	ImportGradObj	PlotMoments
AnisoFilterTensor	EigensysCalc	ImportNii	PlotPhyslog
ApplyCrop	EigenvalCalc	ImportNiiDiff	PlotRespiract
AutoCropData	EigenvecCalc	ImportNiiDix	PlotSegmentMask
BayesianIVIMFit2	EnergyCalc	ImportNiiT1	PlotSegments
BayesianIVIMFit3	EPGSignal	ImportNiiT2	PlotSequence
BlochSeries	EPGT2Fit	ImportPhyslog	PlotSimulation
Bmatrix	ErrorPlot	ImportRespirect	PlotSimulationAngle
BmatrixCalc	ExcludeSlices	ImportVol	PlotSimulationAngleI
BmatrixConv	ExpNoZero	InvertDataset	PlotSimulationHist
BmatrixInv	ExportBmat	IVIMCalc	PlotSimulationVec
BmatrixRot	ExportBval	IVIMCorrectData	PlotSpectrum
BmatrixToggle	ExportBvec	IVIMFunction	Pulses
BullseyePlot	ExportNii	IVIMResiduals	QMRIToolsFuncPrint
BvalRead	ExportVol	JoinSets	QMRIToolsFunctions
CalculateGfactor	ExtractNiiFiles	LapFilter	QMRIToolsPackages
CalculateMoments	FACalc	ListSpherePlot	RadialSample
CalculateWallMap	FConvert	LoadCoilSetup	ReadBrukerDiff
CalibrateEPGT2Fit	FConverti	LoadCoilTarget	ReadBvalue
CardiacCoordinateSystem	FiberDensityMap	LoadFiberTracts	ReadDicom
CardiacSegment	FiberLengths	LogNoZero	ReadDicomDiff
CentralAxes	FileSelect	MADNoZero	ReadDicomDir
ClearTemporaryVariables	FinalGrads	MakeCoilLayout	ReadDicomDirDiff
CoilsSNRCalc	FindCoilPosition	MakeECVBloodMask	ReadGradients
ColorFAPlot	FindCrop	MakeNoisePlots	ReadTransformParamet
CompilebleFunctions	FindMaxDimensions	MakeSliceImages	ReadVoxSize
CompressNiiFiles	FindOrder	MakeWeightMask	RegisterCardiacData
ConcatenateDiffusionData	FindOutliers	Mask	RegisterData
ConditionNumberCalc	FitData	MaskData	RegisterDataSplit
ConvertGrads	FracCorrect	MaskHelix	RegisterDataTransfor
Correct	FullGrad	MeanNoZero	RegisterDataTransfor
CorrectBmatrix	GenerateGradients	MeanRange	RegisterDiffusionDat
CorrectGradients	GenerateGradientsGUI	MeanSignal	RegisterDiffusionDat
CorrectJoinSetMotion	GetGradientScanOrder	MeanStd	RemoveIsoImages
CorrectParMap	GetMaskData	MedCouple	RemoveMaskOverlaps
CreateDiffData	GetMaskMeans	MedianNoZero	RescaleData
CreateHeart	GetPulseProfile	MemoryUsage	RescaleSegmentation
CreateT2Dictionary	GetSliceData	MergeSegmentations	ResidualCalc
CropData	GetSliceNormal	NNLeastSquares	ReverseCrop
CutData	GetSliceNormalDir	NoiseCorrelation	RMSNoZero
Data2DToVector	GetSlicePositions	NoiseCovariance	ROIMask
Data3DToVector	GetSpinSystem	NonLinearEPGFit	SaveImage
DataTransformation	GfactorSimulation	NormalizeData	SegmentMask
DatRead	GradBmatrix	NumberTableForm	SequencePulseAcquire
DatTot	GradientPlot	OverPlusCalc	SequenceSpinEcho
DatTotXLS	GradRead	PadToDimensions	SequenceSteam
DatWrite	GradSeq	ParameterCalc	SequenceTSE
DcmToNii	GridData	ParameterFit	SetupDataStructure
DeNoise	GridData3D	ParameterFit2	ShiftPar
Deriv	HelixAngleCalc	PCADeNoise	SigmaCalc
DevideNoZero	Hist	PCAFitEq	Signal
DictionaryMinSearch	Hist2	PCAFitHist	SimAddPhase
DixonReconstruct	HistogramPar	PhaseAlign	SimAngleParameters
DixonToPercent	HomoginizeData	PlotContour	SimEvolve
DriftCorrect	ImportBmat	PlotCorrection	SimHamiltonian
DTItoolExp	ImportBval	PlotData	SimParameters

Options

AffineDirections	FilterShape	Method	PlotSolution
AnisoFilterSteps	FilterSize	Method	PlotSpace
AnisoKappa	FilterType	Method	PlotStyle
AnisoStepTime	FindTransform	Method	PositiveZ
AnisoWeightType	FitConstraints	Method	PrintTempDirectory
AxesLabel	FitFunction	MethodReg	RadialSamples
AxesMethod	FitOutput	MethodReg	ReadoutBandwidth
BackgroundValue	FitSigma	MethodRegA	ReadoutOutput
BinaryType	FixPseudoDiff	MonitorCalc	ReadoutPhase
BloodMaskRange	FixPseudoDiffSD	MonitorEPGFit	ReadoutSamples
BmatrixOut	FlipAxes	MonitorIVIMCalc	RegistrationTarget
BsplineDirections	FlipBvec	MonitorUnwrap	Reject
BsplineSpacing	FlipGrad	MotionCorrectSets	Reject
BullPlotMethod	FullOutput	NiiDataType	RejectMap
CenterFrequency	FullSphere	NiiMethod	ReportFits
ChainSteps	GetMaskOutput	NiiScaling	Resolutions
CoilArrayPlot	GOutput	NoiseSize	ResolutionsA
CoilSurfaceVoxelSize	GradType	NormalizeIVIM	ReverseData
ColorFunction	GRregularization	NormalizeSets	ReverseSets
ColorFunction	GridLineSpacing	NormalizeSignal	RobustFit
ColorValue	HelixMethod	NumberSamples	RobustFitParameters
CompressNii	HistogramBins	NumberSamplesA	RotateGradient
ConditionCalc	HistogramBinsA	OrderSpan	RotateGradients
ContourStyle	ImageLegend	OutlierIncludeZero	RotationCorrect
ConvertDcm	ImageResolution	OutlierIterations	RowSize
CorrectPar	ImageSize	OutlierMethod	Runs
CropInit	ImageSize	OutlierOutput	SampleStep
CropOutput	ImageSize	OutlierRange	ScaleCorrect
CropPadding	ImageSize	OutputCalibration	Scaling
CutOffMethod	ImageSize	OutputCheckImage	SeedDensity
DeleteTempDirectory	InterpolationOrder	OutputCoilSurface	ShowFit
DeNoiseIterations	InterpolationOrder	OutputImage	ShowHelixPlot
DeNoiseKernel	InterpolationOrderReg	OutputMethod	SliceRange
DeNoiseMonitor	InterpolationOrderRegA	OutputPlot	SliceRangeSamples
DictB1Range	Iterations	OutputSamples	SmartMaskOutput
DictT2fRange	IterationsA	OutputSNR	SmartMethod
DictT2Range	IVIMComponents	OutputTransformation	SmoothHelix
DistanceMeasure	IVIMConstrained	OutputType	SmoothSNR
Distribution	IVIMConstraints	OutputWeights	SortVecs
DixonAmplitudes	IVIMFixed	PadDirection	SpectrumColor
DixonFieldStrength	IVIMTensFit	PaddOverlap	SphereColor
DixonFilterInput	JoinSetSplit	PadValue	SphereSize
DixonFilterOutput	LineStep	Parallelize	SplitMethod
DixonFilterSize	LineThreshold	Parallelize	StartPoints
DixonFrequencies	Linewidth	PCAComponents	StartSlices
DixonIterations	LinewidthShape	PCAFitParameters	Steps
DixonMaskThreshold	MagnetizationVector	PCAKernel	StepSizeI
DixonPrecessions	MakeCheckPlot	PCAOutput	Strictness
DixonTolerance	MaskClosing	PCATolerance	SwitchAxes
DropSamples	MaskCompartment	PCAWeighting	TableAlignments
DropSlices	MaskComponents	PeakNumber	TableDepth
EPGCalibrate	MaskFiltKernel	PhaseEncoding	TableDirections
EPGFitPoints	MaskSmoothing	PlotColor	TableHeadings
EPGMethod	MaskWallMap	PlotLabel	TableMethod
EPGRelaxPars	MeanMethod	PlotLabel	TableSpacing
EPGSmoothB1	MeanRes	PlotRange	TempDirectory
FieldStrength	Method	PlotRange	TensOutput
FileType	Method	PlotRange	TextOffset
FilterMaps	Method	PlotRange	TextSize
			UnitMulti

In[12]:= QMRIToolsFunctions["All", 8]

CardiacTools

Functions

BullseyePlot	HelixAngleCalc
CalculateWallMap	MakeECVBloodMask
CardiacCoordinateSystem	MaskHelix
CardiacSegment	PlotSegmentMask
CentralAxes	PlotSegments
CreateHeart	RadialSample
ECVCalc	TransmuralPlot
ExcludeSlices	

Options

AxesMethod	GridLineSpacing	PlotLabel	StartPoints
BackgroundValue	HelixMethod	PlotRange	StartSlices
BloodMaskRange	ImageSize	PlotStyle	TextOffset
BullPlotMethod	LineStep	RadialSamples	TextSize
ColorFunction	LineThreshold	RowSize	
CutOffMethod	MaskWallMap	ShowFit	
DistanceMeasure	Method	ShowHelixPlot	
DropSamples	OutputCheckImage	SmoothHelix	

CoilTools

Functions

CoilSNRCalc	NoiseCovariance
FindCoilPosition	
LoadCoilSetup	
LoadCoilTarget	
MakeCoilLayout	
MakeNoisePlots	
MakeWeightMask	
NoiseCorrelation	

Options

CoilArrayPlot
CoilSurfaceVoxelSize
ColorFunction
ImageSize
OutputCoilSurface
PlotRange

DenoiseTools

Functions

AnisoFilterTensor
DeNoise
PCADeNoise
PCAFitEq
PCAFitHist
WeightMapCalc

Options

AnisoFilterSteps	Method
AnisoKappa	PCAFitParameters
AnisoStepTime	PCAKernel
AnisoWeightType	PCAOOutput
DeNoiseIterations	PCATolerance
DeNoiseKernel	PCAWeighting
DeNoiseMonitor	PlotSolution
FitSigma	

DixonTools

Functions

DixonReconstruct
DixonToPercent
SimulateDixonSignal
Unwrap
UnwrapSplit

Options

DixonAmplitudes	DixonPrecessions
DixonFieldStrength	DixonTolerance
DixonFilterInput	MonitorUnwrap
DixonFilterOutput	UnwrapDimension
DixonFilterSize	
DixonFrequencies	
DixonIterations	
DixonMaskThreshold	

ElastixTools

Functions

ReadTransformParameters	TransformData
RegisterCardiacData	
RegisterData	
RegisterDataSplit	
RegisterDataTransform	
RegisterDataTransformSplit	
RegisterDiffusionData	
RegisterDiffusionDataSplit	

Options

AffineDirections	InterpolationOrderRegA	OutputTransformation	UseGPU
BsplineDirections	Iterations	PCAComponents	
BsplineSpacing	IterationsA	PrintTempDirectory	
DeleteTempDirectory	MethodReg	RegistrationTarget	
FindTransform	MethodRegA	Resolutions	
HistogramBins	NumberSamples	ResolutionsA	
HistogramBinsA	NumberSamplesA	SplitMethod	
InterpolationOrderReg	OutputImage	TempDirectory	

GeneralTools

Functions

ApplyCrop	DevideNoZero	LogNoZero	QMRIToolsFunctions	SumM
AutoCropData	ExpNoZero	MADNoZero	QMRIToolsPackages	Ten:
ClearTemporaryVariables	FileSelect	MeanNoZero	RescaleData	Ten:
CompilebleFunctions	FindCrop	MedianNoZero	ReverseCrop	Trai
CropData	FindMaxDimensions	MemoryUsage	RMSNoZero	Vec:
CutData	GridData	NNLeastSquares	SaveImage	
Data2DToVector	GridData3D	PadToDimensions	StdFilter	
Data3DToVector	LapFilter	QMRIToolsFuncPrint	StichData	

Options

CropInit	PadDirection
CropOutput	PadValue
CropPadding	WindowTitle
FileType	
ImageResolution	
ImageSize	
InterpolationOrder	
OutputWeights	

GradientTools

Functions

Bmatrix	ConvertGrads	GenerateGradientsGUI	UniqueBvalPosition
BmatrixCalc	CorrectBmatrix	GetGradientScanOrder	
BmatrixConv	CorrectGradients	GetSliceNormal	
BmatrixInv	EnergyCalc	GetSliceNormalDir	
BmatrixRot	FinalGrads	GradBmatrix	
BmatrixToggle	FindOrder	GradSeq	
CalculateMoments	FullGrad	ImportGradObj	
ConditionNumberCalc	GenerateGradients	OverPlusCalc	

Options

ConditionCalc	OutputPlot	UseGrad
FlipAxes	OutputType	VisualOpt
FlipGrad	PhaseEncoding	
FullSphere	Runs	
GradType	Steps	
Method	StepSizeI	
MethodReg	SwitchAxes	
OrderSpan	UnitMulti	

ImportTools

Functions

BvalRead	ReadGradients
GradRead	ReadVoxSize
ReadBrukerDiff	ShiftPar
ReadBvalue	
ReadDicom	
ReadDicomDiff	
ReadDicomDir	
ReadDicomDirDiff	

Options

BmatrixOut
ConvertDcm
RotateGradient
ScaleCorrect

IVIMTools

Functions

BayesianIVIMFit2	IVIMCorrectData
BayesianIVIMFit3	IVIMFunction
CorrectParMap	IVIMResiduals
FConvert	ThetaConv
FConverti	ThetaConvi
FracCorrect	
HistogramPar	
IVIMCalc	

Options

ChainSteps	IVIMComponents	Parallelize
CorrectPar	IVIMConstrained	UpdateStep
FilterMaps	IVIMConstrains	
FilterSize	IVIMFixed	
FilterType	IVIMTensFit	
FitConstrains	Method	
FixPseudoDiff	MonitorIVIMCalc	
FixPseudoDiffSD	OutputSamples	

JcouplingTools

Functions

GetSpinSystem	SimEvolve
PhaseAlign	SimHamiltonian
PlotSpectrum	SimReadout
SequencePulseAcquire	SimRotate
SequenceSpinEcho	SimSignal
SequenceSteam	SimSpoil
SequenceTSE	SysTable
SimAddPhase	

Options

CenterFrequency	ReadoutSamples
FieldStrength	SpectrumColor
Linewidth	
LinewidthShape	
PlotRange	
ReadoutBandwidth	
ReadoutOutput	
ReadoutPhase	

MaskingTools

Functions

GetMaskData	RescaleSegmentation
HomoginizeData	ROIMask
Mask	SegmentMask
MaskData	SmoothMask
MeanSignal	SmoothSegmentation
MergeSegmentations	SplitSegmentations
NormalizeData	
RemoveMaskOverlaps	

Options

GetMaskOutput
MaskClosing
MaskComponents
MaskFiltKernel
MaskSmoothing
UseMask

NiftiTools

Functions

CompressNiiFiles	ImportBval	ImportNiiT2
DcmToNii	ImportBvalvec	
ExportBmat	ImportBvec	
ExportBval	ImportExploreDTItens	
ExportBvec	ImportNii	
ExportNii	ImportNiiDiff	
ExtractNiiFiles	ImportNiiDix	
ImportBmat	ImportNiiT1	

Options

CompressNii
FlipBvec
Method
NiiDataType
NiiMethod
NiiScaling
RotateGradients

PhysiologyTools

Functions

AlignRespLog
ImportPhyslog
ImportRespirect
PlotPhyslog
PlotRespiract

Options

OutputMethod
SampleStep

PlottingTools

Functions

GetSliceData	PlotData3D
GetSlicePositions	PlotDefGrid
GradientPlot	PlotDuty
ListSpherePlot	PlotIVIM
MakeSliceImages	PlotMoments
PlotContour	PlotSequence
PlotCorrection	
PlotData	

Options

ColorFunction	PeakNumber
ContourStyle	PlotColor
DropSlices	PlotRange
ImageLegend	PlotSpace
ImageSize	PositiveZ
MakeCheckPlot	SphereColor
Method	SphereSize
NormalizeIVIM	

ProcessingTools

Functions

CorrectJoinSetMotion	FitData	MedCouple	SplitSets
DataTransformation	GetMaskMeans	NumberTableForm	
DatTot	Hist	ParameterFit	
DatTotXLS	Hist2	ParameterFit2	
ErrorPlot	InvertDataset	SetupDataStructure	
FiberDensityMap	JoinSets	SmartMask	
FiberLengths	MeanRange	SNRCalc	
FindOutliers	MeanStd	SNRMapCalc	

Options

AxesLabel	MeanMethod	OutlierRange	SmartMaskOutput	TableMethod
ColorValue	Method	OutputSNR	SmartMethod	TableSpacing
FitFunction	MotionCorrectSets	PaddOverlap	SmoothSNR	
FitOutput	NormalizeSets	PlotLabel	Strictness	
ImageSize	OutlierIncludeZero	ReverseData	TableAlignments	
InterpolationOrder	OutlierIterations	ReverseSets	TableDepth	
JoinSetSplit	OutlierMethod	Scaling	TableDirections	
MaskCompartment	OutlierOutput	SeedDensity	TableHeadings	

RelaxometryTools

Functions

CalibrateEPGT2Fit	T2Fit
CreateT2Dictionary	TriExponentialT2Fit
DictionaryMinSearch	
EPGSignal	
EPGT2Fit	
NonLinearEPGFit	
T1Fit	
T1rhoFit	

Options

DictB1Range	Method
DictT2fRange	MonitorEPGFit
DictT2Range	OutputCalibration
EPGCalibrate	
EPGFitPoints	
EPGMethod	
EPGRelaxPars	
EPGSmoothB1	

SimulationTools

Functions

AddNoise	PlotSimulationAngleHist	Tensor
BlochSeries	PlotSimulationHist	
CalculateGfactor	PlotSimulationVec	
CreateDiffData	Pulses	
GetPulseProfile	Signal	
GfactorSimulation	SimAngleParameters	
PlotSimulation	SimParameters	
PlotSimulationAngle	SimulateSliceEPG	

Options

GOutput	SliceRangeSamples
GRegularization	SortVecs
MagnetizationVector	TensOutput
NoiseSize	
PlotRange	
Reject	
ReportFits	
SliceRange	

TensorTools

Functions

ADCCalc	ECalc	SigmaCalc
AngleCalc	EigensysCalc	SortDiffusionData
AngleMap	EigenvalCalc	TensorCalc
ColorFAPlot	EigenvecCalc	TensorCorrect
ConcatenateDiffusionData	FACalc	
Correct	ParameterCalc	
Deriv	RemoveIsoImages	
DriftCorrect	ResidualCalc	

Options

Distribution	Reject
FilterShape	RejectMap
FullOutput	RobustFit
MeanRes	RobustFitParameters
Method	RotationCorrect
MonitorCalc	UseMask
NormalizeSignal	
Parallelize	

VisteTools

Functions

DatRead	ImportVol
DatWrite	LoadFiberTracts
DTItoolExp	
DTItoolExpFile	
DTItoolExpInd	
DTItoolExpTens	
ExportVol	
ImportDTI	

Options

BinaryType

```
In[7]:= QMRIToolsFuncPrint[]
```

CardiacTools

Functions

BullseyePlot[data, segmask] generates a AHA–17 segment bullseye plot.

BullseyePlot[list] generates a AHA–17 segment
bullseye plot of the lists (which needs to have 17 values) provide.

data is a 3D volume used for the plot.

segmask is the AHA–17 segmentation

resulting from the CardiacSegment function when AHA17 is selected.

Output is a bullseye plot or a plotwindow, depending

on the Method which can be "Dynamic" else it will be static. >>

```
Attributes[BullseyePlot] = {Protected, ReadProtected}
```

```
Options[BullseyePlot] = {TextOffset → 0.5, TextSize → 12, PlotRange → Automatic,  
ColorFunction → TemperatureMap, BullPlotMethod → Dynamic, ImageSize → 200}
```

```
SyntaxInformation[BullseyePlot] = {ArgumentsPattern → {_, _., OptionsPattern[]}}
```

CalculateWallMap[mask,vox] calculates the wall distance map and the wall derivative.

Output is {wallmap, wallDerivative}. >>

```
Attributes[CalculateWallMap] = {Protected, ReadProtected}
```

```
Options[CalculateWallMap] = {ShowFit → True, MaskWallMap → True}
```

```
SyntaxInformation[CalculateWallMap] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

QMRITools`CardiacTools`CardiacCoordinateSystem

```
Attributes[CardiacCoordinateSystem] = {Protected, ReadProtected}
```

```
Options[CardiacCoordinateSystem] = {ShowFit → False}
```

```
SyntaxInformation[CardiacCoordinateSystem] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

CardiacSegment[data, mask, off] allows to segment the heart in 1,
4, 6 or AHA-17 segments for each slice 360 radial samples are generated.

data is a background image on which all overlays are projected.

mask is the mask of the left ventricle (same as

used for CentralAxes) and defines the area in which the data is sampled.

off is the centerpoints generated by CentralAxes.

Output is {segmask, segang, {points, slices}}. >>

```
Attributes[CardiacSegment] = {Protected, ReadProtected}
```

```
Options[CardiacSegment] =  
{StartPoints → Default, StartSlices → Default, LineThreshold → 0.25, LineStep → 0.5}
```

```
SyntaxInformation[CardiacSegment] = {ArgumentsPattern → {_, _, _, OptionsPattern[]}}
```

CentralAxes[mask, vox] calculates the center of the

lumen from a mask of the left ventricle. vox is the voxels size, {slice, x, y}.

CentralAxes[mask, maskp, vox] allows for fancy visualization of the other structures using maskp.

Output is {centerpoints, normalvecs, inout} or {centerpoints, normalvecs, inout, fit}. >>

```
Attributes[CentralAxes] = {Protected, ReadProtected}
```

```
Options[CentralAxes] = {ShowFit → True, RowSize → Automatic, AxesMethod → Cubic}
```

```
SyntaxInformation[CentralAxes] = {ArgumentsPattern → {_, _, _., OptionsPattern[]}}
```

CreateHeart[] creates a simulated left ventricle shape.

CreateHeart[pars] creates a simulated left ventricle shape with predefined parameters pars.

Output is the heart shape, the voxel size and

the parameters needed to generate the heart, {mask, vox, pars}.

```
Attributes[CreateHeart] = {Protected, ReadProtected}
```

```
SyntaxInformation[CreateHeart] = {ArgumentsPattern → {_.}}
```

ECVCalc[T1pre, T1post, hema] calculates the ECVmap using MakeECVBloodMask.

ECVCalc[T1pre, T1post, bloodMask, hema] calculates the ECVmap using bloodMask. >>

```
Attributes[ECVCalc] = {Protected, ReadProtected}
```

ExcludeSlices[data] excludes slices that

do not look like the others based on various distance measures.

Output is an array with 1 or 0 with the dimensions {slices, diff dirs} >>

```
Attributes[ExcludeSlices] = {Protected, ReadProtected}
```

```
Options[ExcludeSlices] = {CutOffMethod → Auto, DistanceMeasure → 5}
```

```
SyntaxInformation[ExcludeSlices] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

HelixAngleCalc[eigenvectors, mask, vox] calculates

the helix angle matrix of cardiac data using only a left ventricle mask.

HelixAngleCalc[eigenvectors, mask, vox] calculates the helix angle matrix of cardiac data using only a left ventricle mask, and a maskp for visualization.

HelixAngleCalc[eigenvectors, mask, centerpoint, vec, inout, vox] calculates the helix angle matrix of cardiac data using only a left ventricle mask.

HelixAngleCalc[eigenvectors, mask, maskp, centerpoint, vec, inout, vox] calculates the helix angle matrix of cardiac data using a left ventricle mask and a maskp for visualization.

eigenvectors are the tensor eigenvectors calculated with EigenvecCalc.

mask is a mask of the left ventricle.

maskp is a mask used for visualization.

vox is the voxels size, {slice, x, y}.

The following values are calculated

automaticlay Using CentralAxes but can also be provided as an input.

centerpoint is the center of each slice calculated with CentralAxes.

inout is the inner and outer radius calculated with CentralAxes.

vec is the vector describing the central axes of the heart, calculated with CentralAxes.

Output is the fiber angle matrix FAM = {9, slice, x, y} or {FAM, plot}. >>

```
Attributes[HelixAngleCalc] = {Protected, ReadProtected}
```

```
Options[HelixAngleCalc] = {ShowHelixPlot → True, HelixMethod → Slow, AxesMethod → Quadratic}
```

```
SyntaxInformation[HelixAngleCalc] = {ArgumentsPattern → {_, _, _, _., OptionsPattern[]}}
```

MakeECVBloodMask[T1pre, T1post] makes a bloodpool mask based on the T1pre and T1post images. It assumes that the heart is cropped with the blood in the center. >>

```
Attributes[MakeECVBloodMask] = {Protected, ReadProtected}
```

```
Options[MakeECVBloodMask] = {BloodMaskRange → {1400, {0, 700}}, OutputCheckImage → True}
```

```
SyntaxInformation[MakeECVBloodMask] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

MaskHelix[helix, mask] masks helix angle data, sets

the background to -100 and allows for Median filter of the helix mask.

helix can be a single map or the FAM.

Output is the masked helix angle data. >>

```
Attributes[MaskHelix] = {Protected, ReadProtected}
```

```
Options[MaskHelix] = {BackgroundValue → -100, SmoothHelix → False}
```

```
SyntaxInformation[MaskHelix] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

`PlotSegmentMask[mask, segmask, vox]` plots the mask segments created by `CardiacSegment`.

`mask` is a mask the left ventricle that was used in the `CardiacSegment`.

`segmask` is the output of `CardiacSegment`.

`vox` is the voxels size, {slice, x, y}.

Output is a plot window. >>

```
Attributes[PlotSegmentMask] = {Protected, ReadProtected}
```

```
SyntaxInformation[PlotSegmentMask] = {ArgumentsPattern -> {_, _, _}}
```

`PlotSegments[mask, data, segang]` shows how the heart will be sampled by `RadialSample`.

`mask` is a mask the left ventricle that was used in the `CardiacSegment`.

function and the `segang` is the output of the `CardiacSegmentFunction`.

Output is a plot window. >>

```
Attributes[PlotSegments] = {Protected, ReadProtected}
```

```
Options[PlotSegments] = {RadialSamples -> 10}
```

```
SyntaxInformation[PlotSegments] = {ArgumentsPattern -> {_, _, _, OptionsPattern[]}}
```

`RadialSample[mask, data, segang]` radially samples the provided parametermap data.

The `mask` should be a mask of the left ventricle that was used in the `CardiacSegment`.

`segang` is the output of the `CardiacSegmentFunction`.

Output is {points, vals} which are ordered as indicated by the user. >>

```
Attributes[RadialSample] = {Protected, ReadProtected}
```

```
Options[RadialSample] = {RadialSamples -> 10, DropSamples -> 0}
```

```
SyntaxInformation[RadialSample] = {ArgumentsPattern -> {_, _, _, OptionsPattern[]}}
```

`TransmuralPlot[data]` plots transmural profiles of the data which are created by `RadialSample`.

`data` can be a single profile or a list of profiles.

In the second case the mean and standard deviations are plotted.

Output is a plot of the transmural profile. >>


```
Attributes[TransmuralPlot] = {Protected, ReadProtected}
```

```
Options[TransmuralPlot] = {GridLineSpacing → 10, PlotStyle → RGBColor[1, 0, 0],  
  PlotRange → Automatic, ImageSize → 300, Method → Median, PlotLabel → None}
```

```
SyntaxInformation[TransmuralPlot] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

Options

AxesMethod is an option for HelixAngleCalc and CentralAxes. Can be "Linear", "Quadratic", "Cubic". >>

```
Attributes[AxesMethod] = {Protected, ReadProtected}
```

BackgroundValue is an option for MaskHelix. Sets the background value (default is -100). >>

```
Attributes[BackgroundValue] = {Protected, ReadProtected}
```

BloodMaskRange is an option for MakeECVBloodMask. >>

```
Attributes[BloodMaskRange] = {Protected, ReadProtected}
```

BullPlotMethod is an option for BullseyePlot. Can be "Dynamic" or "Normal".
"Dynamic" allows to change plotting parameters in Manipulation window. >>

```
Attributes[BullPlotMethod] = {Protected, ReadProtected}
```

ColorFunction is an option for graphics functions
that specifies a function to apply to determine colors of elements. >>

```
Attributes[ColorFunction] = {Protected}
```

CutOffMethod is an option for ExcludeSlices. Default
value is "Auto" or it can be a fixed percentage (value between 0 and .5) >>

```
Attributes[CutOffMethod] = {Protected, ReadProtected}
```

DistanceMeasure is an option for ExcludeSlices. Default value is 5. (1 ManhattanDistance, 2
SquaredEuclideanDistance, 3 EuclideanDistance, 4 Correlation, 5 SpearmanRho >>

```
Attributes[DistanceMeasure] = {Protected, ReadProtected}
```

DropSamples is an option for RadialSample and PlotSegments. Defines how many samples
are dropped from star and end. Can be a number or set (start, end) of numbers. >>

```
Attributes[DropSamples] = {Protected, ReadProtected}
```

GridLineSpacing is an option of TransmuralPlot. It defines the spacing of the gridlines. >>

```
Attributes[GridLineSpacing] = {Protected, ReadProtected}
```

HelixMethod is an option for HelixAngleCalc. Can be "Slow" or "Fast".
 "Slow" uses wall distance interpolation and can take long for high res datasets.
 "Fast" uses wall distance calculation using circular approximation of the ventricle. >>

```
Attributes[HelixMethod] = {Protected, ReadProtected}
```

ImageSize is an option that specifies the overall size of an image to display for an object. >>

```
Attributes[ImageSize] = {Protected}
```

LineStep is an option for CardiacSegment. >>

```
Attributes[LineStep] = {Protected, ReadProtected}
```

LineThreshold is an option for CardiacSegment. Can be number between
 0 and 1. Increasing the value will decrease the amount of wall sampled. >>

```
Attributes[LineThreshold] = {Protected, ReadProtected}
```

MaskWallMap is an option for CalculateWallMap. if True or False. >>

```
Attributes[MaskWallMap] = {Protected, ReadProtected}
```

Method is an option for various algorithm-intensive
 functions that specifies what internal methods they should use. >>

```
Attributes[Method] = {Protected}
```

QMRITools`CardiacTools`OutputCheckImage

```
Attributes[OutputCheckImage] = {Protected, ReadProtected}
```

PlotLabel is an option for graphics functions that specifies an overall label for a plot. >>

```
Attributes[PlotLabel] = {Protected}
```

PlotRange is an option for graphics functions
 that specifies what range of coordinates to include in a plot. >>

```
Attributes[PlotRange] = {Protected, ReadProtected}
```

PlotStyle is an option for plotting and related
 functions that specifies styles in which objects are to be drawn. >>

```
Attributes[PlotStyle] = {Protected}
```

RadialSamples is an option for RadialSample and
 PlotSegments. Defines how many transmural samples are taken. >>

```
Attributes[RadialSamples] = {Protected, ReadProtected}
```

RowSize is an option for CentralAxes. defines the number of images per showing the segmentation.
Can be "Automatic" or an integer. >>

Attributes[RowSize] = {Protected, ReadProtected}

ShowFit is an option for CentralAxes. True shows the fit of the central axes. >>

Attributes[ShowFit] = {Protected, ReadProtected}

ShowHelixPlot is an option for HelixAngleCalc. If true the it
also outputs a visualization of the local myocardial coordinate system. >>

Attributes[ShowHelixPlot] = {Protected, ReadProtected}

SmoothHelix is an option for MaskHelix, sets the kernelsize for the MedianFilter. >>

Attributes[SmoothHelix] = {Protected, ReadProtected}

StartPoints is an option for CardiacSegment.
Value is "Default" or the point list given by CardiacSegment. >>

Attributes[StartPoints] = {Protected, ReadProtected}

StartSlices is an option for CardiacSegment. Value is "Default" or the list given by CardiacSegment. >>

Attributes[StartSlices] = {Protected, ReadProtected}

TextOffset is an option for BullseyePlot. Determines where the text is placed, can be 0 to 1. >>

Attributes[TextOffset] = {Protected, ReadProtected}

TextSize is an option for BullseyePlot. Determines the text size. >>

Attributes[TextSize] = {Protected, ReadProtected}

CoilTools

Functions

CoilSNRCalc[coils, noise] calculates the sensitivity
weighted snr of multiple coil elements using magnitude signal and noise.

Output is {data, noise, sos, snr, sigmap, weights}. >>

Attributes[CoilSNRCalc] = {Protected, ReadProtected}

SyntaxInformation[CoilSNRCalc] = {ArgumentsPattern → {_, _}}

`FindCoilPosition[weights]` finds the coil position by locating the highest intensity location in the coil weight map, which can be obtained by `LoadCoilSetup` or `SumOfSquares`. Internally it uses `MakeWeightMask` to remove the noise of the weightmaps. `FindCoilPosition[weights, mask]` limits the search region to the provided mask. >>

```
Attributes[FindCoilPosition] = {Protected, ReadProtected}
```

```
Options[FindCoilPosition] = {OutputCoilSurface → False, CoilSurfaceVoxelSize → {1, 1, 1}}
```

```
SyntaxInformation[FindCoilPosition] = {ArgumentsPattern → {_, _., OptionsPattern[]}}
```

`LoadCoilSetup[file]` load a very specific type of coil experiment, a dynamic scan with a setup of which the second dynamic is a noise measurement. The input file is the Nii file that contains the individually reconstructed coil images and the noise data. Internally it uses `CoilSNRCalc` and `SumOfSquares`.

Output is the coil data with coil noise data and snrmap based on the `SumOfSquares` addition, the SOS reconstruction and the SOS weights. `{dataC, noiseC, sosC, snrC, sigmapC, weights, vox}`. >>

```
Attributes[LoadCoilSetup] = {Protected, ReadProtected}
```

```
SyntaxInformation[LoadCoilSetup] = {ArgumentsPattern → {_}}
```

`LoadCoilTarget[file]` loads a very specific type of experiment, a dynamic scan with the second dynamic is a noise measurement. The input file is the Nii file that contains the scanner reconstruction and the noise data. Internally it uses `SNRMapCalc`,

Output is the reconstructed data with noise data and snrMap `{dataC, noiseC, sosC, snrC, sigmapC, weights, vox}`. >>

```
Attributes[LoadCoilTarget] = {Protected, ReadProtected}
```

```
SyntaxInformation[LoadCoilTarget] = {ArgumentsPattern → {_}}
```

`MakeCoilLayout[{name, size, number}]` makes a coil grid with label name, partitioned in size rows and with label number. `MakeCoilLayout[{name, size, number}, val]` makes a coil grid with label name, partitioned in size rows and with label the val at location number. `MakeCoilLayout[{coils..}]` same but for multiple coils grids. Each coil grid is defined as `{name, size, number}`. `MakeCoilLayout[{coils..}, val]` same but for multiple coil grids. >>

```
Attributes[MakeCoilLayout] = {Protected, ReadProtected}
```

```
Options[MakeCoilLayout] = {PlotRange → Automatic, ColorFunction → SunsetColors, ImageSize → 100, CoilArrayPlot → False}
```

MakeNoisePlots[noise] returns a grid of plots of the noise per channel
 MakeNoisePlots[noise, {met, prt}] met can be "Grid" with
 prt a number or Automatic. Else all plots will be returned as a list of plots.
 MakeNoisePlots[noise, {met, prt}, sub] sub defines how much the noise is
 subsampled, default is 40 (every 40th sample is used in plot). >>

```
Attributes[MakeNoisePlots] = {Protected, ReadProtected}
```

```
SyntaxInformation[MakeNoisePlots] = {ArgumentsPattern -> {_, _., OptionsPattern[]}}
```

MakeWeightMask[weights] creates a mask
 of homogeneous regions of weightmaps removing the noise. >>

```
Attributes[MakeWeightMask] = {Protected, ReadProtected}
```

```
SyntaxInformation[MakeWeightMask] = {ArgumentsPattern -> {_}}
```

NoiseCorrelation[noise] calculates the noise correlation matrix, noise is {nrCoils, noise Samples}. >>

```
Attributes[NoiseCorrelation] = {Protected, ReadProtected}
```

```
SyntaxInformation[NoiseCorrelation] = {ArgumentsPattern -> {_}}
```

NoiseCovariance[noise] calculates the noise covariance matrix, noise is {nrCoils, noise Samples}. >>

```
Attributes[NoiseCovariance] = {Protected, ReadProtected}
```

```
SyntaxInformation[NoiseCovariance] = {ArgumentsPattern -> {_}}
```

Options

CoilArrayPlot is an option for MakeCoilLayout. If True
 and values are provided it makes an arrayplot of the coil layouts >>

```
Attributes[CoilArrayPlot] = {Protected, ReadProtected}
```

CoilSurfaceVoxelSize is an option for
 FindCoilPosition. Specifies the voxel size used for OutputCoilSurface. >>

```
Attributes[CoilSurfaceVoxelSize] = {Protected, ReadProtected}
```

ColorFunction is an option for graphics functions
 that specifies a function to apply to determine colors of elements. >>

```
Attributes[ColorFunction] = {Protected}
```

ImageSize is an option that specifies the overall size of an image to display for an object. >>

```
Attributes[ImageSize] = {Protected}
```

OutputCoilSurface is an option for FindCoilPosition. If
 set true it will also output a SurfacePlot of the coil location volume. >>

```
Attributes[OutputCoilSurface] = {Protected, ReadProtected}
```

PlotRange is an option for graphics functions
that specifies what range of coordinates to include in a plot. >>

```
Attributes[PlotRange] = {Protected, ReadProtected}
```

DenoiseTools

Functions

AnisoFilterTensor[tens, diffdata] Filter the
tensor tens using an anisotropic diffusion filter (Perona–Malik).
It uses the diffusion weighted data diffdata to find edges that are not visible in the tensor.
Edge weights based on the diffusion data are averaged over all normalized diffusion direction.

Output is the smoothed tensor. >>

```
Attributes[AnisoFilterTensor] = {Protected, ReadProtected}
```

```
Options[AnisoFilterTensor] =  
{AnisoWeightType → 2, AnisoKappa → 5., AnisoStepTime → 1, AnisoFilterSteps → 5}
```

```
SyntaxInformation[AnisoFilterTensor] = {ArgumentsPattern → {_, _}, OptionsPattern[] }
```

DeNoise[data, sigma, filtersize] removes Rician noise with standard deviation "sigma"
from the given dataset using a kernel with size "filtersize" a gaussian kernel.
DeNoise[data, sigma, filtersize, Kernel → "kerneltype"] removes Rician noise with standard deviation
"sigma" from the given dataset using a kernel with size "filtersize" and type "kerneltype".

Output is data denoised. >>

```
Attributes[DeNoise] = {Protected, ReadProtected}
```

```
Options[DeNoise] = {DeNoiseKernel → Gaussian, DeNoiseMonitor → False, DeNoiseIterations → 1}
```

```
SyntaxInformation[DeNoise] = {ArgumentsPattern → {_, _, _}, OptionsPattern[] }
```

PCADeNoise[data] removes rician noise from the data with PCA.
PCADeNoise[data, mask] removes rician noise from the data with PCA only withing the mask.
PCADeNoise[data, mask, sig] removes rician noise from the data
with PCA only withing the mask using sig as prior knowledge or fixed value.

Output is de {data denoise, sigma map} by default if PCAOutput is Full
then fitted {data dnoise, {sigma fit, average sigma}, {number components,
number of fitted voxesl, number of max fits}, total fit –time per 500 ittt}. >>

```
Attributes[PCADeNoise] = {Protected, ReadProtected}
```

```
Options[PCADeNoise] = {PCAKernel → 5, PCAFitParameters → {10, 6, 10}, FitSigma → False,  
  PCAOutput → Full, Method → Equation, PCATolerance → 0, PCAWeighting → True}
```

```
SyntaxInformation[PCADeNoise] = {ArgumentsPattern → {_, __, __, OptionsPattern[]}}
```

PCAFitEq[data] fits the marchencopasteur distribution to the PCA of the data using grid search.
PCAFitEq[data, sig] fits the marchencopasteur distribution to
the PCA of the data using sig as start value or fixed value using grid search.

Output is {simga, number of noise comp, and denoised matrix}. >>

```
Attributes[PCAFitEq] = {Protected, ReadProtected}
```

```
SyntaxInformation[PCAFitEq] = {ArgumentsPattern → {_, __, __, __}}
```

PCAFitHist[data] fits the marchencopasteur distribution to the PCA of the data using hist fit.
PCAFitHist[data, sig] fits the marchencopasteur distribution
to the PCA of the data using sig as start value or fixed value using hist fit.

Output is {simga, number of noise comp, and denoised matrix, itterations}. >>

```
Attributes[PCAFitHist] = {Protected, ReadProtected}
```

```
Options[PCAFitHist] = {PlotSolution → False, FitSigma → True, PCAFitParameters → {10, 6, 10}}
```

```
SyntaxInformation[PCAFitHist] = {ArgumentsPattern → {_, __, OptionsPattern[]}}
```

WeightMapCalc[diffdata] calculates a weight map which is used in AnisoFilterTensor.

Output is a weight map of the diffdata which is high in isotropic regions and low at edges. >>

```
Attributes[WeightMapCalc] = {Protected, ReadProtected}
```

```
Options[WeightMapCalc] = {AnisoWeightType → 2, AnisoKappa → 10.}
```

```
SyntaxInformation[WeightMapCalc] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

Options

AnisoFilterSteps is an option for AnisoFilterTensor and
defines the amoutn of diffusin steps taken. Higher is more smoothing >>

```
Attributes[AnisoFilterSteps] = {Protected, ReadProtected}
```

AnisoKappa is an option for AnisoFilterTensor and WeightMapCalc and
defines the weighting strenght, all data is normalize to 100 before filetering. >>

```
Attributes[AnisoKappa] = {Protected, ReadProtected}
```

AnisoStepTime is an option for AnisoFilterTensor
and defines the diffusion time, when small more step are needed. >>

Attributes[AnisoStepTime] = {Protected, ReadProtected}

AnisoWeightType is an option for AnisoFilterTensor and WeightMapCalc and defines
the weighting, either 1, the exponent of $(-g/\kappa)$ or 2, $1/(1+g/\kappa)$. >>

Attributes[AnisoWeightType] = {Protected, ReadProtected}

DeNoiseIterations is and option for DeNoise. Specifies the number of the denoising iterations. >>

Attributes[DeNoiseIterations] = {Protected, ReadProtected}

DeNoiseKernel is and option for DeNoise. Values can be "Disk", "Box" or "Gaussian". >>

Attributes[DeNoiseKernel] = {Protected, ReadProtected}

DeNoiseMonitor is and option for DeNoise. Monitor the denoising progres. >>

Attributes[DeNoiseMonitor] = {Protected, ReadProtected}

FitSigma is an option of PCAFitHist, PCAFitEq and
PCADeNoise, if set True sig is fitted if set False sigma is fixed to input value. >>

Attributes[FitSigma] = {Protected, ReadProtected}

Method is an option for various algorithm-intensive
functions that specifies what internal methods they should use. >>

Attributes[Method] = {Protected}

PCAFitParameters is an option of PCADeNoise and PCAFitHist. {nb, pi,
maxit} = bins, initial signal components, maximum number of itterations. >>

Attributes[PCAFitParameters] = {Protected, ReadProtected}

PCAKernel is an option of PCADeNoise. It sets the kernel size. >>

Attributes[PCAKernel] = {Protected, ReadProtected}

PCAOOutput is an option of PCADeNoise. If output is full the
output is {datao, {output[[1]], sigmat}, {output[[2]], output[[3]], j}, timetot}.
Else the output is {datao, sigmat}. >>

Attributes[PCAOOutput] = {Protected, ReadProtected}

PCATolerance is an option of PCADeNoise and shuld be an integer > 0.
Default value is 0. When increased the denoise method removes less noise. >>

Attributes[PCATolerance] = {Protected, ReadProtected}

PCAWeighting is an option of PCADeNoise and can be True of False. Default value is False. When True the weights of the per voxel result are calculated based on the number of non noise components. >>

Attributes[PCAWeighting] = {Protected, ReadProtected}

PlotSolution is an option for PCAFitHist, if set true it displays the fitting iterations. >>

Attributes[PlotSolution] = {Protected, ReadProtected}

DixonTools

Functions

DixonReconstruct[real, imag, echo] reconstructs Dixon data with initial guess $B_0 = 0$ and $T_2^* = 0$.

DixonReconstruct[real, imag, echo, B_0] reconstructs Dixon data with initial guess $T_2^* = 0$.

DixonReconstruct[real, imag, echo, B_0 , t2] reconstructs Dixon data.

real is the real data in radians.

imag is the imaginary data in radians.

B_0 can be estimated from two phase images using Unwrap.

T2 can be estimated from multiple echos using T2fit.

Output is {{watF,fatF},{watSig,fatSig},{inphase,outphase},{ B_0 , T_2^* },iterations}. >>

Attributes[DixonReconstruct] = {Protected, ReadProtected}

Options[DixonReconstruct] = {DixonPrecessions → -1, DixonFieldStrength → 3,
DixonFrequencies → {{0}, {3.8, 3.4, 3.13, 2.67, 2.46, 1.92, 0.57, -0.6}},
DixonAmplitudes → {{1}, {0.089, 0.598, 0.047, 0.077, 0.052, 0.011, 0.035, 0.066}},
DixonIterations → 5, DixonTolerance → 0.01, DixonMaskThreshold → 0.05,
DixonFilterInput → True, DixonFilterOutput → True, DixonFilterSize → 1}

SyntaxInformation[DixonReconstruct] = {ArgumentsPattern → {_, _, _, _., _., OptionsPattern[]}}

DixonToPercent[water, fat] converts the dixon water and fat data to percent maps.

Output is {waterFraction, fatFraction}. >>

Attributes[DixonToPercent] = {Protected, ReadProtected}

SyntaxInformation[DixonToPercent] = {ArgumentsPattern → {_, _}}

SimulateDixonSignal[echo, fr, B_0 , T2] simulates an Dixon gradient echo sequence with echotimes echo in ms, fat fraction fr, field of resonance B_0 in Hz and relaxation T2 in ms.

```
Attributes[SimulateDixonSignal] = {Protected, ReadProtected}
```

```
Options[SimulateDixonSignal] = {DixonPrecessions → -1, DixonFieldStrength → 3,  
  DixonFrequencies → {{0}, {3.8, 3.4, 3.13, 2.67, 2.46, 1.92, 0.57, -0.6}},  
  DixonAmplitudes → {{1}, {0.089, 0.598, 0.047, 0.077, 0.052, 0.011, 0.035, 0.066}}}
```

```
SyntaxInformation[SimulateDixonSignal] = {ArgumentsPattern → {_, _, _, _}, OptionsPattern[]}
```

Unwrap[data] unwraps the given dataset. >>

```
Attributes[Unwrap] = {Protected, ReadProtected}
```

```
Options[Unwrap] = {MonitorUnwrap → True, UnwrapDimension → 2D}
```

```
SyntaxInformation[Unwrap] = {ArgumentsPattern → {_, _}, OptionsPattern[]}
```

UnwrapSplit[phase, data] unwraps the give phase dataset but splits the data into left and right using SplitData based in the data and performs the unwrapping seperately. >>

```
Attributes[UnwrapSplit] = {Protected, ReadProtected}
```

```
Options[UnwrapSplit] = {}
```

```
SyntaxInformation[UnwrapSplit] = {ArgumentsPattern → {_, _}, OptionsPattern[]}
```

Options

DixonAmplitudes is an options for DixonReconstruct. Defines the amplitudes of the fat peaks being used. >>

```
Attributes[DixonAmplitudes] = {Protected, ReadProtected}
```

DixonFieldStrength is an options for DixonReconstruct. Defines the fieldstrengths on which the data was acquired. >>

```
Attributes[DixonFieldStrength] = {Protected, ReadProtected}
```

DixonFilterInput is an options for DixonReconstruct. If True the input b0 and T2star values are smoothed using a gaussian kernel. >>

```
Attributes[DixonFilterInput] = {Protected, ReadProtected}
```

DixonFilterOutput is an options for DixonReconstruct. If True the out b0 and T2star values are smoothed Median filter and lowpassfiltering after which the water and fat maps are recomputed. >>

```
Attributes[DixonFilterOutput] = {Protected, ReadProtected}
```

DixonFilterSize is an options for DixonReconstruct. Defines the number of voxel with which the input b0 and T2star values are smoothed. >>

```
Attributes[DixonFilterSize] = {Protected, ReadProtected}
```

DixonFrequencies is an options for
DixonReconstruct. Defines the frequencies of the fat peaks being used. >>

Attributes[DixonFrequencies] = {Protected, ReadProtected}

DixonIterations is an options for DixonReconstruct. Defines the maximum iterations the fit can use. >>

Attributes[DixonIterations] = {Protected, ReadProtected}

DixonMaskThreshold is an options for DixonReconstruct. Defines at which threshold the dixon reconstruction considers a voxel to be background noise. Default values is 0.05. >>

Attributes[DixonMaskThreshold] = {Protected, ReadProtected}

DixonPrecessions is an options for
DixonReconstruct. Defines the rotation of the signal {-1,1} default is -1. >>

Attributes[DixonPrecessions] = {Protected, ReadProtected}

DixonTolerance is an options for DixonReconstruct. Defines at which change per iteration of b0 and R2star the iterative methods stops. Default value is 0.1. >>

Attributes[DixonTolerance] = {Protected, ReadProtected}

MonitorUnwrap is an option for Unwrap. Monitor the unwrapping progress. >>

Attributes[MonitorUnwrap] = {Protected, ReadProtected}

UnwrapDimension is an option for Unwrap. Can be "2D" or "3D". 2D is for unwrapping 2D images or unwrapping the individual images from a 3D dataset (does not unwrap in the slice direction). 3D unwraps a 3D dataset in all dimensions. >>

Attributes[UnwrapDimension] = {Protected, ReadProtected}

ElastixTools

Functions

ReadTransformParameters[directory] reads the tranformation parameters generated by RegisterData. The directory should be the TempDirectory were the registration is stored. DeleteTempDirectory should be False.

Output is the affine transformation vector per volume. >>

Attributes[ReadTransformParameters] = {Protected, ReadProtected}

SyntaxInformation[ReadTransformParameters] = {ArgumentsPattern -> {_}}

RegisterCardiacData[data] registers the data using a 2D algorithm. data can be 3D or 4D.
 RegisterCardiacData[{data,vox}] registers the data series using the given voxel size.
 RegisterCardiacData[{data,mask}] registers the data series only using data within the mask.
 RegisterCardiacData[{data,mask,vox}] registers the
 data series using the given voxel size only using data within the mask.

Output is the registered data. >>

```
Attributes[RegisterCardiacData] = {Protected, ReadProtected}
```

```
Options[RegisterCardiacData] = {RegistrationTarget -> Mean,  
  Iterations -> 250, Resolutions -> 1, HistogramBins -> 64, NumberSamples -> 2000,  
  InterpolationOrderReg -> 3, BsplineSpacing -> 30, BsplineDirections -> {1, 1, 1},  
  AffineDirections -> {1, 1, 1}, MethodReg -> affine, OutputImage -> True,  
  TempDirectory -> Default, DeleteTempDirectory -> True, PrintTempDirectory -> True,  
  OutputTransformation -> False, UseGPU -> {False, Automatic}, PCAComponents -> 3}
```

```
SyntaxInformation[RegisterCardiacData] = {ArgumentsPattern -> {_, OptionsPattern[]}}
```

RegisterData[data] registers the data series. If data is 3D it performs multiple 2D registration, if data is 4D it performs multiple 3D registration. The target is the first image or volume in the series.

RegisterData[{data, vox}] registers the data series using the given voxel size.

RegisterData[{data, mask}] registers the data series only using data within the mask.

RegisterData[{data, mask, vox}] registers the data series using the given voxel size only using data within the mask.

RegisterData[target, moving] registers the moving data to the target data. target can be 2D or 3D. moving can be the same of one dimension higher than the target.

RegisterData[{target, mask, vox}, {moving, mask, vox}] registers the data using the given voxel size only using data within the mask.

RegisterData[{target, vox}, moving] registers the data using the given voxel size.

RegisterData[target, {moving, vox}] registers the data using the given voxel size.

RegisterData[{target, vox}, {moving, vox}] registers the data using the given voxel size.

RegisterData[{target, mask}, moving] registers the data series only using data within the mask.

RegisterData[target, {moving, mask}] registers the data series only using data within the mask.

RegisterData[{target, mask}, moving] registers the data series only using data within the mask.

RegisterData[{target, mask}, {moving, mask}] registers the data series only using data within the mask.

RegisterData[target, {moving, mask, vox}] registers the data series using the given voxel size only using data within the mask.

RegisterData[{target, mask}, {moving, mask, vox}] registers the data series using the given voxel size only using data within the mask.

RegisterData[{target, vox}, {moving, mask, vox}] registers the data series using the given voxel size only using data within the mask.

RegisterData[{target, mask, vox}, moving] registers the data series using the given voxel size only using data within the mask.

RegisterData[{target, mask, vox}, {moving, mask}] registers the data series using the given voxel size only using data within the mask.

RegisterData[{target, mask, vox}, {moving, vox}] registers the data series using the given voxel size only using data within the mask.

RegisterData[{target, mask}, {moving, vox}] registers the data series using the given voxel size only using data within the mask.

RegisterData[{target, vox}, {moving, mask}] registers the data series using the given voxel size only using data within the mask.

Output is the registered data with the dimensions of the moving data.

If OutputTransformation is True it also outputs the translation, rotation scale and skew of all images or volumes. >>

```
Attributes[RegisterData] = {Protected, ReadProtected}
```

```
Options[RegisterData] =
```

```
{Iterations → 250, Resolutions → 1, HistogramBins → 64, NumberSamples → 2000,
 InterpolationOrderReg → 3, BsplineSpacing → 30, BsplineDirections → {1, 1, 1},
 AffineDirections → {1, 1, 1}, MethodReg → affine, OutputImage → True,
 TempDirectory → Default, DeleteTempDirectory → True, PrintTempDirectory → True,
 OutputTransformation → False, UseGPU → {False, Automatic}, PCAComponents → 3}
```

```
SyntaxInformation[RegisterData] = {ArgumentsPattern → {_, _}, OptionsPattern[] }
```

RegisterDataSplit[target, moving] is identical to RegisterData
data however left and right side of the data are registered separately.

Splitting the data is done using the function CutData and merged with Stich data.

Output is the registered data. >>

```
Attributes[RegisterDataSplit] = {Protected, ReadProtected}
```

```
Options[RegisterDataSplit] =
```

```
{Iterations → 250, Resolutions → 1, HistogramBins → 64, NumberSamples → 2000,
 InterpolationOrderReg → 3, BsplineSpacing → 30, BsplineDirections → {1, 1, 1},
 AffineDirections → {1, 1, 1}, MethodReg → affine, OutputImage → True, TempDirectory → Default,
 DeleteTempDirectory → True, PrintTempDirectory → True, OutputTransformation → False,
 UseGPU → {False, Automatic}, PCAComponents → 3, SplitMethod → Mean}
```

```
SyntaxInformation[RegisterDataSplit] = {ArgumentsPattern → {_, _}, OptionsPattern[] }
```

RegisterDataTransform[target, moving, {moving2nd,
vox}] performs the registration exactly as RegisterData. target and
moving are the inputs for Registerdata, which can be {data,mask,vox}.

After the registration is done the moving2nd data is deformed
according to the output of the registration of moving.

moving2nd can have the same dimensions
of moving or one dimension higher (e.g. 3D and 3D or 3D and 4D).

Output is {registered moving, deformed moving2nd}. >>

```
Attributes[RegisterDataTransform] = {Protected, ReadProtected}
```

```
Options[RegisterDataTransform] =
```

```
{Iterations → 250, Resolutions → 1, HistogramBins → 64, NumberSamples → 2000,
 InterpolationOrderReg → 3, BsplineSpacing → 30, BsplineDirections → {1, 1, 1},
 AffineDirections → {1, 1, 1}, MethodReg → affine, OutputImage → True,
 TempDirectory → Default, DeleteTempDirectory → True, PrintTempDirectory → True,
 OutputTransformation → False, UseGPU → {False, Automatic}, PCAComponents → 3}
```

```
SyntaxInformation[RegisterDataTransform] = {ArgumentsPattern → {_, _, _}, OptionsPattern[] }
```

RegisterDataTransformSplit[target, moving, {moving2nd, vox}] is identical to RegisterDataTransform with the same functionality as RegisterDataSplit. This means the data is split in two using the function CutData and merged with Stich data.

Output is {registered moving, deformed moving2nd}. >>

```
Attributes[RegisterDataTransformSplit] = {Protected, ReadProtected}
```

```
Options[RegisterDataTransformSplit] =
{Iterations → 250, Resolutions → 1, HistogramBins → 64, NumberSamples → 2000,
 InterpolationOrderReg → 3, BsplineSpacing → 30, BsplineDirections → {1, 1, 1},
 AffineDirections → {1, 1, 1}, MethodReg → affine, OutputImage → True, TempDirectory → Default,
 DeleteTempDirectory → True, PrintTempDirectory → True, OutputTransformation → False,
 UseGPU → {False, Automatic}, PCAComponents → 3, SplitMethod → Mean}
```

```
SyntaxInformation[RegisterDataTransformSplit] =
{ArgumentsPattern → {_, _, _}, OptionsPattern[]}
```

RegisterDiffusionData[{dtidata, vox}] registers a diffusion dataset.

dtidata should be 4D {slice, diff, x, y}. vox is the voxel size of the data.

RegisterDiffusionData[{dtidata, dtimask, vox}] registers the data series using the given voxel size only using data within the mask.

RegisterDiffusionData[{dtidata, vox}, {anatdata, vox}] registers a diffusion dataset. The diffusion data is also registered to the anatdata.

RegisterDiffusionData[{dtidata, dtimask, vox}, {anatdata, vox}] registers the data series using the given voxel size only using data within the mask.

RegisterDiffusionData[{dtidata, vox}, {anatdata, anatmask, vox}] registers the data series using the given voxel size only using data within the mask.

RegisterDiffusionData[{dtidata, dtimask, vox}, {anatdata, anatmask, vox}] registers the data series using the given voxel size only using data within the mask.

Output is the registered dtidata and, if anatdata is given, the registered dtidata in anatomical space. If OutputTransformation is True it also outputs the translation, rotation scale and skew of all images or volumes. >>

```
Attributes[RegisterDiffusionData] = {Protected, ReadProtected}
```

```
Options[RegisterDiffusionData] =
{Iterations → 250, Resolutions → 1, HistogramBins → 64, NumberSamples → 2000,
 InterpolationOrderReg → 3, BsplineSpacing → 30, BsplineDirections → {0, 1, 1},
 AffineDirections → {0, 1, 1}, MethodReg → affineDTI, OutputImage → True,
 TempDirectory → Default, DeleteTempDirectory → True, PrintTempDirectory → True,
 OutputTransformation → False, UseGPU → {False, Automatic}, PCAComponents → 3,
 IterationsA → 1000, ResolutionsA → 1, HistogramBinsA → 64, NumberSamplesA → 20000,
 InterpolationOrderRegA → 1, MethodRegA → {rigidDTI, bspline}, RegistrationTarget → F1st}
```

```
SyntaxInformation[RegisterDiffusionData] = {ArgumentsPattern → {_, _}, OptionsPattern[]}
```

RegisterDiffusionDataSplit[dtidata, vox] is identical to Register diffusion data however left and right side of the data are registered separately.

RegisterDiffusionDataSplit[{dtidata, vox}, {anatdata, vox}] is identical to Register diffusion data however left and right side of the data are registered separately.

RegisterDiffusionDataSplit[{dtidata, dtimask, vox}, {anatdata, anatmask, vox}] is identical to Register diffusion data however left and right side of the data are registered separately.

Splitting the data is done using the function CutData and merged with Stich data.

Output is the registered data. >>

```
Attributes[RegisterDiffusionDataSplit] = {Protected, ReadProtected}
```

```
Options[RegisterDiffusionDataSplit] := Options[RegisterDiffusionData]
```

```
SyntaxInformation[RegisterDiffusionDataSplit] = {ArgumentsPattern -> {_, _., OptionsPattern[]}}
```

TransformData[{data, vox}] deforms the data according to the last output of register data.

The directory should be the TempDirectory where the registration is stored. DeleteTempDirectory should be False. >>

```
Attributes[TransformData] = {Protected, ReadProtected}
```

```
Options[TransformData] = {TempDirectory -> Default,
  FindTransform -> Auto, DeleteTempDirectory -> All, PrintTempDirectory -> True}
```

```
SyntaxInformation[TransformData] = {ArgumentsPattern -> {_, OptionsPattern[]}}
```

Options

AffineDirections is an option for RegisterData and RegisterDiffusionData.

It gives the directions in which data can be moved when registering diffusion data to anatomical space. >>

```
Attributes[AffineDirections] = {Protected, ReadProtected}
```

BsplineDirections is an option for RegisterData and RegisterDiffusionData.

It gives the direction in which the bsplines are allowed to move when registering diffusion data to anatomical space. >>

```
Attributes[BsplineDirections] = {Protected, ReadProtected}
```

BsplineSpacing is an option for RegisterData, RegisterDiffusionData, RegisterCardiacData and RegisterDataTransform.

It specifies the spacing of the bsplines if the method is "bspline". >>

```
Attributes[BsplineSpacing] = {Protected, ReadProtected}
```


DeleteTempDirectory an options for RegisterData, RegisterDiffusionData, RegisterCardiacData and RegisterDataTransform. It specifies if the temp directory should be deleted after the registration is finished. >>

```
Attributes[DeleteTempDirectory] = {Protected, ReadProtected}
```

FindTransform is an option for TransformData and RegisterTransformData. It specifies where to find the transformfile. >>

```
Attributes[FindTransform] = {Protected, ReadProtected}
```

HistogramBins is an options for RegisterData, RegisterDiffusionData, and RegisterDataTransform. It specifies the number of bins of the joined histogram used by the registration functions. >>

```
Attributes[HistogramBins] = {Protected, ReadProtected}
```

HistogramBinsA is an option for RegisterDiffusionData. It specifies the number of bins of the joined histogram used when registering diffusion data to anatomical space. >>

```
Attributes[HistogramBinsA] = {Protected, ReadProtected}
```

InterpolationOrderReg is an options for RegisterData, RegisterDiffusionData, and RegisterDataTransform. It specifies the interpolation order used in the registration functions. >>

```
Attributes[InterpolationOrderReg] = {Protected, ReadProtected}
```

InterpolationOrderRegA is an option for RegisterDiffusionData. It specifies the interpolation order used in the registration functions when registering diffusion data to anatomical space. >>

```
Attributes[InterpolationOrderRegA] = {Protected, ReadProtected}
```

Iterations is an options for RegisterData, RegisterDiffusionData, and RegisterDataTransform. It specifies the number of iterations used by the registration functions. >>

```
Attributes[Iterations] = {Protected, ReadProtected}
```

IterationsA is an option for RegisterDiffusionData. It specifies the number of iterations used when registering diffusion data to anatomical space. >>

```
Attributes[IterationsA] = {Protected, ReadProtected}
```

MethodReg is an options for RegisterData, RegisterDiffusionData, RegisterCardiacData and RegisterDataTransform. It specifies which registration method to use. Mehtods can be be "rigid", "affine", "bspline" or "cyclyc". >>

```
Attributes[MethodReg] = {Protected, ReadProtected}
```

MethodRegA is an option for RegisterDiffusionData. It specifies which registration method to use when registering diffusion data to anatomical space. Mehtods can be be "rigid", "affine" or "bspline". >>

```
Attributes[MethodRegA] = {Protected, ReadProtected}
```

NumberSamples is an options for RegisterData, RegisterDiffusionData, and RegisterDataTransform. It specifies the number of random samples that are taken each iteration used by the registration functions. >>

```
Attributes[NumberSamples] = {Protected, ReadProtected}
```

NumberSamplesA is an option for RegisterDiffusionData. It specifies the number of random samples that are taken each iteration when registering diffusion data to anatomical space. >>

```
Attributes[NumberSamplesA] = {Protected, ReadProtected}
```

OutputImage is an options for RegisterData, RegisterDiffusionData, and RegisterDataTransform. It specifies if the result image should be written in the TempDirectory as nii file. >>

```
Attributes[OutputImage] = {Protected, ReadProtected}
```

OutputTransformation is an option for RegisterData ad RegisterDiffusionData. It specifies if the tranformation paramters (translation, rotation, scale and skew) should be given as output in the registration functions. >>

```
Attributes[OutputTransformation] = {Protected, ReadProtected}
```

PCAComponents is an option for RegisterData. It speciefies how many PCA components are used if method is set to "PCA" >>

```
Attributes[PCAComponents] = {Protected, ReadProtected}
```

PrintTempDirectory is an options for RegisterData, RegisterDiffusionData, RegisterCardiacData and RegisterDataTransform. It specifies if the location of the temp directory should be deployed. >>

```
Attributes[PrintTempDirectory] = {Protected, ReadProtected}
```

RegistrationTarget is an option for RegisterDiffusionData and RegisterCardiacData. Specifies which target to use for registration if using "rigid", "affine" or "bspline" as MethodReg. If the MethodReg is "cyclic" or "PCA" it does not need a target and this option does nothing. Values can be "First", "Mean" or "Median". >>

```
Attributes[RegistrationTarget] = {Protected, ReadProtected}
```

Resolutions is an option for RegisterData, RegisterDiffusionData, and RegisterDataTransform. It specifies the number of scale space resolutions used by the registration functions. >>

```
Attributes[Resolutions] = {Protected, ReadProtected}
```

ResolutionsA is an option for RegisterDiffusionData. It specifies the number of scale space resolutions used when registering diffusion data to anatomical space. >>

```
Attributes[ResolutionsA] = {Protected, ReadProtected}
```

SplitMethod is an option for RegisterDataSplit and RegisterDataTransformSplit. values can be "mean", "moving", "target" >>

```
Attributes[SplitMethod] = {Protected, ReadProtected}
```

TempDirectory is an option for RegisterData, RegisterDiffusionData, RegisterCardiacData and RegisterDataTransform. It specifies the temporary directory used to perform and output the registration. >>

```
Attributes[TempDirectory] = {Protected, ReadProtected}
```

UseGPU is an option for RegisterData. The value is {bool, gpu} where bool is True or False, and gpu is the gpu ID which is an integer or Automatic. >>

```
Attributes[UseGPU] = {Protected, ReadProtected}
```

GeneralTools

Functions

ApplyCrop[data,crop] applies the cropped region obtained from CropData to the data.
 ApplyCrop[data,crop,{voxorig,voxnew}]
 applies the cropped region obtained from CropData to the data. >>

```
Attributes[ApplyCrop] = {Protected, ReadProtected}
```

```
SyntaxInformation[ApplyCrop] = {ArgumentsPattern -> {_, _, _..}}
```

AutoCropData[data] crops the data by removing all background zeros.
 AutoCropData[data,pad] crops the data by removing all background zeros with padding of pad. >>

```
Attributes[AutoCropData] = {Protected, ReadProtected}
```

```
Options[AutoCropData] = {CropPadding → 5}
```

```
SyntaxInformation[AutoCropData] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

ClearTemporaryVariables[] Clear temporary variables. >>

```
Attributes[ClearTemporaryVariables] = {Protected, ReadProtected}
```

```
SyntaxInformation[ClearTemporaryVariables] = {ArgumentsPattern → {_.}}
```

CompilebleFunctions[] generates a list of all compilable functions. >>

```
Attributes[CompilebleFunctions] = {Protected, ReadProtected}
```

```
SyntaxInformation[CompilebleFunctions] = {ArgumentsPattern → {}}
```

CropData[data] creates a dialog window to crop the data (assumes voxsize (1,1,1)).

CropData[data,vox] creates a dialog window to crop the data. >>

```
Attributes[CropData] = {Protected, ReadProtected}
```

```
Options[CropData] = {CropOutput → All, CropInit → Automatic}
```

```
SyntaxInformation[CropData] = {ArgumentsPattern → {_, _., OptionsPattern[]}}
```

CutData[data] splits the data in two equal sets left and right. >>

```
Attributes[CutData] = {Protected, ReadProtected}
```

```
SyntaxInformation[CutData] = {ArgumentsPattern → {_, _.,}}
```

Data2DToVector[data] convert the data to vector.

Data2DToVector[data,mask] convert the data within the mask to vector.

the data can be reconstructed using VectorToData.

output is the vectorized data and a list containing the original data dimensions and a list with the data coordinates. {vec, {dim,pos}}. >>

```
Attributes[Data2DToVector] = {Protected, ReadProtected}
```

```
SyntaxInformation[Data2DToVector] = {ArgumentsPattern → {_, _.,}}
```

Data3DToVector[data] convert the data to vector..

Data3DToVector[data,mask] convert the data within the mask to vector.

the data can be reconstructed using VectorToData.

output is the vectorized data and a list containing the original data dimensions and a list with the data coordinates. {vec, {dim,pos}}. >>

```
Attributes[Data3DToVector] = {Protected, ReadProtected}
```

```
SyntaxInformation[Data3DToVector] = {ArgumentsPattern → {_, _}}
```

DevideNoZero[a, b] devides a/b but when b=0 the result is 0. a can be a number or vector. >>

```
Attributes[DevideNoZero] = {Protected, ReadProtected}
```

```
SyntaxInformation[DevideNoZero] = {ArgumentsPattern → {_, _}}
```

ExpNoZero[val] return the Exp of the val which can be anny dimonsion array. if val=0 the output is 0. >>

```
Attributes[ExpNoZero] = {Protected, ReadProtected}
```

```
SyntaxInformation[ExpNoZero] = {ArgumentsPattern → {_}}
```

FileSelect[action] creates a systemdialog wicht returs

file/foldername action can be "FileOpen", "FileSave" or "Directory".

FileSelect[action, {type}] same but allows the definition of filetypes

for "FileOpen" and "FileSave" e.g. "jpg" or "pdf". >>

```
Attributes[FileSelect] = {Protected, ReadProtected}
```

```
Options[FileSelect] = {WindowTitle → Automatic}
```

```
SyntaxInformation[FileSelect] = {ArgumentsPattern → {_, __, __, OptionsPattern[]}}
```

FindCrop[data] finds the crop values of the data by removing all zeros surrounding the data. >>

```
Attributes[FindCrop] = {Protected, ReadProtected}
```

```
Options[FindCrop] = {CropPadding → 5}
```

```
SyntaxInformation[FindCrop] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

FindMaxDimensions[{data1, data2, ...}] finds

the maximal dimensions of all datasets. Each dataset is 3D. >>

```
Attributes[FindMaxDimensions] = {Protected, ReadProtected}
```

```
SyntaxInformation[FindMaxDimensions] = {ArgumentsPattern → {_}}
```

GridData[{data1, data2, ...}, part] makes a grid of multiple datasets with part sets on each row >>

```
Attributes[GridData] = {Protected, ReadProtected}
```

```
SyntaxInformation[GridData] = {ArgumentsPattern → {_, _}}
```

GridData3D[{data1, data2, ...}, part] same as grid data, but only

works on 4D data where the data is gridded in axial, coronal and sagital.

```
Attributes[GridData3D] = {Protected, ReadProtected}
```

LapFilter[data] Laplacian filter of data with kernel size 0.8.

LapFilter[data, ker] Laplacian filter of data with kernel ker. >>

Attributes[LapFilter] = {Protected, ReadProtected}

LogNoZero[val] return the log of the val which can be any dimension array. if val=0 the output is 0. >>

Attributes[LogNoZero] = {Protected, ReadProtected}

SyntaxInformation[LogNoZero] = {ArgumentsPattern → {_}}

MADNoZero[vec] return the MAD error of the vec which can

be any dimension array. if vec={0...} the output is 0. Zeros are ignored >>

Attributes[MADNoZero] = {Protected, ReadProtected}

SyntaxInformation[MADNoZero] = {ArgumentsPattern → {_}}

MeanNoZero[data] calculates the mean of the data ignoring the zeros. >>

Attributes[MeanNoZero] = {Protected, ReadProtected}

SyntaxInformation[MeanNoZero] = {ArgumentsPattern → {_, _}}

MedianNoZero[data] calculates the Median of the data ignoring the zeros. >>

Attributes[MedianNoZero] = {Protected, ReadProtected}

SyntaxInformation[MedianNoZero] = {ArgumentsPattern → {_, _}}

MemoryUsage[] gives a table of which definitions use up memory.

MemoryUsage[n] gives a table of which

definitions use up memory, where n is the amount of definitions to show. >>

Attributes[MemoryUsage] = {Protected, ReadProtected}

SyntaxInformation[MemoryUsage] = {ArgumentsPattern → {_.}}

NNLeastSquares[A, y] performs a Non Negative Linear Least Squares fit.

finds an x that solves the linear least-squares problem for the matrix equation $Ax=y$.

output is the solution x. >>

Attributes[NNLeastSquares] = {Protected, ReadProtected}

SyntaxInformation[NNLeastSquares] = {ArgumentsPattern → {_, _}}

PadToDimensions[data, dim] pads the data to dimensions dim. >>

Attributes[PadToDimensions] = {Protected, ReadProtected}

Options[PadToDimensions] = {PadValue → 0., PadDirection → Center}

SyntaxInformation[PadToDimensions] = {ArgumentsPattern → {_, _, OptionsPattern[]}}

QMRIToolsFuncPrint[] gives a list of all the QMRITools functions with their usage information. >>

```
Attributes[QMRIToolsFuncPrint] = {Protected, ReadProtected}
```

```
SyntaxInformation[QMRIToolsFuncPrint] = {ArgumentsPattern -> {_.}}
```

QMRIToolsFunctions[] give list of all the QMRITools packages, functions and options.

QMRIToolsFunctions[p] print a table with length p of all the QMRITools functions and options.

QMRIToolsFunctions["toolbox"] gives a list of all the functions and options in toolbox.

QMRIToolsFunctions["toolbox", p] gives a table of length p of all the functions and options in toolbox. If toolbox is "All" it will list all toolboxes. >>

```
Attributes[QMRIToolsFunctions] = {Protected, ReadProtected}
```

```
SyntaxInformation[QMRIToolsFunctions] = {ArgumentsPattern -> {_, _}}
```

QMRIToolsPackages[] give list of all the QMRITools packages. >>

```
Attributes[QMRIToolsPackages] = {Protected, ReadProtected}
```

```
SyntaxInformation[QMRIToolsPackages] = {ArgumentsPattern -> {}}
```

RescaleData[data,dim] rescales image/data to given dimensions.

RescaleData[data,{vox1, vox2}] rescales image/data from size vox1 to size vox2. >>

```
Attributes[RescaleData] = {Protected, ReadProtected}
```

```
Options[RescaleData] = {InterpolationOrder -> 3}
```

```
SyntaxInformation[RescaleData] = {ArgumentsPattern -> {_, _, OptionsPattern[]}}
```

ReverseCrop[data,dim,crop] reverses the crop

on the cropped data with crop values crop to the original size dim.

ReverseCrop[data,dim,crop,{voorig,voxnew}] reverses the crop on the cropped data with crop values crop to the original size dim. >>

```
Attributes[ReverseCrop] = {Protected, ReadProtected}
```

```
SyntaxInformation[ReverseCrop] = {ArgumentsPattern -> {_, _, _, _}}
```

RMSNoZero[vec] return the RMS error of the vec which can

be any dimension array. if vec={0...} the output is 0. Zeros are ignored >>

```
Attributes[RMSNoZero] = {Protected, ReadProtected}
```

```
SyntaxInformation[RMSNoZero] = {ArgumentsPattern -> {_}}
```

SaveImage[image] exports graph to image,

ImageSize, FileType and ImageResolution can be given as options.

SaveImage[image, "filename"] exports graph to image with "filename",

ImageSize, FileType and ImageResolution can be given as options. >>

```
Attributes[SaveImage] = {Protected, ReadProtected}
```

```
Options[SaveImage] = {ImageSize → 6000, FileType → .jpg, ImageResolution → 300}
```

```
SyntaxInformation[SaveImage] = {ArgumentsPattern → {_, _}, OptionsPattern[] }
```

StdFilter[data] StandardDeviation filter of data using gaussian kernel 2.

StdFilter[data, ker] StandardDeviation filter of data using kernel with size ker. >>

```
Attributes[StdFilter] = {Protected, ReadProtected}
```

StichData[datall,datarr] joins left and right part of the data generated by CutData. >>

```
Attributes[StichData] = {Protected, ReadProtected}
```

```
SyntaxInformation[StichData] = {ArgumentsPattern → {_, _} }
```

SumOfSquares[{data1, data2, ..., datan}] calculates the sum of squares of the datasets.

Output is the SoS and the weights, or just the SoS. >>

```
Attributes[SumOfSquares] = {Protected, ReadProtected}
```

```
Options[SumOfSquares] = {OutputWeights → True}
```

```
SyntaxInformation[SumOfSquares] = {ArgumentsPattern → {_, OptionsPattern[] } }
```

TensMat[tensor] transforms tensor form vector format

{xx,yy,zz,xy,xz,yz} to matrix format {{xx,xy,xz},{xy,yy,yz},{xz,yz,zz}}. >>

```
Attributes[TensMat] = {Protected, ReadProtected}
```

```
SyntaxInformation[TensMat] = {ArgumentsPattern → {_} }
```

TensVec[tensor] transforms tensor form matrix format

{{xx,xy,xz},{xy,yy,yz},{xz,yz,zz}} to vector format {xx,yy,zz,xy,xz,yz}. >>

```
Attributes[TensVec] = {Protected, ReadProtected}
```

```
SyntaxInformation[TensVec] = {ArgumentsPattern → {_} }
```

TransData[data,dir] Rotates the dimesions of the

data to left or righthg. For example {z,x,y} to {x,y,z} dir is "l" or "r". >>

```
Attributes[TransData] = {Protected, ReadProtected}
```

```
SyntaxInformation[TransData] = {ArgumentsPattern → {_, _} }
```

VectorToData[vec, {dim,pos}] converts the vectroized data, using

Data2DToVector or Data3DToVector, back to its original Dimensoins >>

```
Attributes[VectorToData] = {Protected, ReadProtected}
```

```
SyntaxInformation[VectorToData] = {ArgumentsPattern → {_, {_, _} } }
```

Options

CropInit is an option for CropData. By default the crop is not initialized but can be with `{{xmin,xmax},{ymin,ymax},{zmin,zmax}}`. >>

```
Attributes[CropInit] = {Protected, ReadProtected}
```

CropOutput is an option for CropData, can be "All", "Data" or "Crop". >>

```
Attributes[CropOutput] = {Protected, ReadProtected}
```

CropPadding is an option for AutoCropData or FindCrop. It specifies how much padding to use around the data. >>

```
Attributes[CropPadding] = {Protected, ReadProtected}
```

FileType["file"] gives the type of a file, typically File, Directory, or None. >>

```
Attributes[FileType] = {Protected}
```

ImageResolution is an option for Export, Rasterize, and related functions that specifies at what resolution bitmap images should be rendered. >>

```
Attributes[ImageResolution] = {Protected}
```

ImageSize is an option that specifies the overall size of an image to display for an object. >>

```
Attributes[ImageSize] = {Protected}
```

InterpolationOrder is an option for Interpolation, as well as ListLinePlot, ListPlot3D, ListContourPlot, and related functions, that specifies what order of interpolation to use. >>

```
Attributes[InterpolationOrder] = {Protected}
```

OutputWeights is an option for SumOfSquares. If True it also outputs the SoS weights. >>

```
Attributes[OutputWeights] = {Protected, ReadProtected}
```

PadDirection is an option for PadToDimensions. It specifies the direction of padding, "Center", "Left" or "Right". >>

```
Attributes[PadDirection] = {Protected, ReadProtected}
```

PadValue is an option for PadToDimensions. It specifies the value of the padding. >>

```
Attributes[PadValue] = {Protected, ReadProtected}
```

WindowTitle is an option that specifies the title to give for a window. >>

```
Attributes[WindowTitle] = {Protected}
```

GradientTools

Functions

Bmatrix[bvec,grad] creates bmatrix from grad and bvec in form {-bxx, -byy, -bzz, -bxy, -bxz, -byz, 1}.
 Bmatrix[{bvec,grad}] creates bmatrix from grad and bvec in form {bxx, byy, bzz, bxy, bxz, byz}. >>

Attributes[Bmatrix] = {Protected, ReadProtected}

Options[Bmatrix] = {Method → DTI}

SyntaxInformation[Bmatrix] = {ArgumentsPattern → {_, _., OptionsPattern[]}}

BmatrixCalc["folder", grads] calculates the true bmatrix from the exported sequence parameters from the philips scanner that are stored in "folder" for each of the gradient directions grads. >>

Attributes[BmatrixCalc] = {Protected, ReadProtected}

Options[BmatrixCalc] = {UseGrad → {1, 1, {1, 1}, 1, 1}, OutputType → Matrix,
 Method → Numerical, StepSizeI → 0.05, UnitMulti → $\frac{1}{1000}$, PhaseEncoding → A,
 FlipAxes → {{1, 1, 1}, {1, 1, 1}}, SwitchAxes → {{1, 2, 3}, {1, 2, 3}}}

SyntaxInformation[BmatrixCalc] = {ArgumentsPattern → {_, _, _., OptionsPattern[]}}

BmatrixConv[bm] converts the bmatrix from 7 to 6 or from 6 to 7. >>

Attributes[BmatrixConv] = {Protected, ReadProtected}

SyntaxInformation[BmatrixConv] = {ArgumentsPattern → {_}}

BmatrixInv[bm] generates a bvecotr and gradiens directions form a given bmatrx.
 BmatrixInv[bm, bvi] generates a bvecotr and
 gradiens directions form a given bmatrx using the given bvalues bvi. >>

Attributes[BmatrixInv] = {Protected, ReadProtected}

SyntaxInformation[BmatrixInv] = {ArgumentsPattern → {_, _}}

BmatrixRot[bmat, rotmat] Rotates the B-matrix. >>

Attributes[BmatrixRot] = {Protected, ReadProtected}

SyntaxInformation[BmatrixRot] = {ArgumentsPattern → {_, _}}

BmatrixToggle[bmat, axes, flip], axes can be any order
 of {"x", "y", "z"}. flip should be {1,1,1},{1,1,-1},{1,-1,1} or {-1,1,1}. >>

Attributes[BmatrixToggle] = {Protected, ReadProtected}

SyntaxInformation[BmatrixToggle] = {ArgumentsPattern → {_, _, _}}

CalculateMoments[{Gt, hw, te}, t] calculates the 0th to 3th order moments
 of the sequence created by GradSeq. Output is {{Gt, M0, M1, M2, M3}, vals}. >>

Attributes[CalculateMoments] = {Protected, ReadProtected}

SyntaxInformation[CalculateMoments] = {ArgumentsPattern → {_, _}}

ConditionNumberCalc[grads] calculates the condition number of the gradient set. >>

```
Attributes[ConditionNumberCalc] = {Protected, ReadProtected}
```

```
SyntaxInformation[ConditionNumberCalc] = {ArgumentsPattern -> { _ }}
```

ConvertGrads[grad, bv] converts the gradients to txt format, which is needed for FinalGrads. >>

```
Attributes[ConvertGrads] = {Protected, ReadProtected}
```

```
SyntaxInformation[ConvertGrads] = {ArgumentsPattern -> { _, _ }}
```

CorrectBmatrix[bmat, transformation] corrects the bmatrix bmat with the tranformation parameters from RegisterData or RegisterDiffusionData.

Output is the corrected bmatrix. >>

```
Attributes[CorrectBmatrix] = {Protected, ReadProtected}
```

```
Options[CorrectBmatrix] = {MethodReg -> Full}
```

```
SyntaxInformation[CorrectBmatrix] = {ArgumentsPattern -> { _, _, OptionsPattern[] }}
```

CorrectGradients[grad, transformation] corrects the gradient directions grad with the tranformation parameters from RegisterData or RegisterDiffusionData.

Output is the corrected gradient vector. >>

```
Attributes[CorrectGradients] = {Protected, ReadProtected}
```

```
Options[CorrectGradients] = {MethodReg -> Rotation}
```

```
SyntaxInformation[CorrectGradients] = {ArgumentsPattern -> { _, _, OptionsPattern[] }}
```

EnergyCalc[grads] calculates the total Energy of the gradient set. >>

```
Attributes[EnergyCalc] = {Protected, ReadProtected}
```

```
SyntaxInformation[EnergyCalc] = {ArgumentsPattern -> { _ }}
```

FinalGrads[grtxt,{int,intn},{rand,order}] finalizes the gradient txt file.
grtxt is the output from the function ConvertGrads, which convert the grad to txt format.
int is True or False, if set to True it interleaves b=0 gradients every intn directions.
rand indicates if the gradients need to be randomized, for this it uses the order which is the output of FindOrder. >>

```
Attributes[FinalGrads] = {Protected, ReadProtected}
```

```
SyntaxInformation[FinalGrads] = {ArgumentsPattern -> { _, { _, _ }, { _, _ } }}
```

FindOrder[grad,bv] finds the optimal order of the gradient directions which minimizes the duty cycle.

The output is needed for FinalGrads.

grad is a list of gradient sets and bv is a

list of b-values with the same number as the list of gradient sets. >>

```
Attributes[FindOrder] = {Protected, ReadProtected}
```

```
Options[FindOrder] = {OrderSpan → Auto}
```

```
SyntaxInformation[FindOrder] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

FullGrad is an option for Grad. Default is True. When true

the gradient directions will be loaded with the first gradient {0,0,0}. >>

```
Attributes[FullGrad] = {Protected, ReadProtected}
```

GenerateGradients[numb] optimizes a set with numb gradients, numb must be an integer.

GenerateGradients[{numb, fixed}] optimizes a set with numb gradients, numb must be an integer and fixed a list of 3D coordinates e.g. {{0,0,1},{0,1,0}}. The fixed gradients will not be moved.

GenerateGradients[{numb1, numb2 ...}, alpha] optimizes a multi shell gradient set with numb gradients per shell. If alpha is set to 0.5 equal importance is given to the optimal distribution of each shell in the entire set. If alpha is 0 only the sub shells will be optimized, if alpha is set to 1 only the global set will be optimized. >>

```
Attributes[GenerateGradients] = {Protected, ReadProtected}
```

```
Options[GenerateGradients] = {Steps → 1000, Runs → 1,  
VisualOpt → False, GradType → Normal, ConditionCalc → False, FullSphere → False}
```

```
SyntaxInformation[GenerateGradients] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

GenerateGradientsGUI[] runs the

GenerateGradients function in GUI with output for the Philips system. >>

```
Attributes[GenerateGradientsGUI] = {Protected, ReadProtected}
```

GetGradientScanOrder[grad, bval] determines

the scanorder based on the txt file provided to the scanner as input.

GetGradientScanOrder[file, grad, bval] determines the scanorder based on the txt file provided to the scanner as input. >>

```
Attributes[GetGradientScanOrder] = {Protected, ReadProtected}
```

```
SyntaxInformation[GetGradientScanOrder] = {ArgumentsPattern → {_, _, _}}
```

GetSliceNormal[file] imports the slice normal from a DICOM image. >>

```
Attributes[GetSliceNormal] = {Protected, ReadProtected}
```

```
SyntaxInformation[GetSliceNormal] = {ArgumentsPattern → {_, _}}
```

GetSliceNormalDir[file] imports the slice normal from a enhanced dicom image. >>

```
Attributes[GetSliceNormalDir] = {Protected, ReadProtected}
```

```
SyntaxInformation[GetSliceNormalDir] = {ArgumentsPattern -> { _ }}
```

GradBmatrix[Gt, hw, te, t] Calculates the true bmatrix from the sequence created by GradSeq. >>

```
Attributes[GradBmatrix] = {Protected, ReadProtected}
```

```
Options[GradBmatrix] = {OutputPlot -> False, Method -> Analytical, StepSizeI -> 0.025}
```

```
SyntaxInformation[GradBmatrix] = {ArgumentsPattern -> { _, _, _, _, OptionsPattern[] }}
```

GradSeq[pars, t, grad] Creates a sequence from the gradient pars imported by ImportGradObj. >>

```
Attributes[GradSeq] = {Protected, ReadProtected}
```

```
Options[GradSeq] = {UseGrad -> {0, 1, {1, 0}, 1}, FlipGrad -> False, UnitMulti -> 1,
  PhaseEncoding -> A, FlipAxes -> {{1, 1, 1}, {1, 1, 1}}, SwitchAxes -> {{1, 2, 3}, {1, 2, 3}}}
```

```
SyntaxInformation[GradSeq] = {ArgumentsPattern -> { _, _, _, OptionsPattern[] }}
```

ImportGradObj[folder] Imports the gradient par files exported from the philips scanner. >>

```
Attributes[ImportGradObj] = {Protected, ReadProtected}
```

```
SyntaxInformation[ImportGradObj] = {ArgumentsPattern -> { _ }}
```

OverPlusCalc[grads] determines the minimal overplus factor of of the gradient set. >>

```
Attributes[OverPlusCalc] = {Protected, ReadProtected}
```

```
SyntaxInformation[OverPlusCalc] = {ArgumentsPattern -> { _ }}
```

UniqueBvalPosition[bval] generates a list of all the unique bvalues and their positions.

UniqueBvalPosition[bval, num] generates a list of all the unique bvalues and their positions that are present in the dataset equal or more than num times >>

```
Attributes[UniqueBvalPosition] = {Protected, ReadProtected}
```

```
SyntaxInformation[UniqueBvalPosition] = {ArgumentsPattern -> { _, _ }}
```

Options

ConditionCalc is an option for GenerateGradients if set to true

GenerateGradients will also give the condition number evolution of the system. >>

```
Attributes[ConditionCalc] = {Protected, ReadProtected}
```

FlipAxes is an option for GradSeq. Defaul value is {{1,1,1},{1,1,1}}. First three

values are for diffusion gradients last three are for the acquisition gradients. >>

```
Attributes[FlipAxes] = {Protected, ReadProtected}
```

FlipGrad is an option for GradSeq. When FlipGrad is true the gr180 is flipped. >>

```
Attributes[FlipGrad] = {Protected, ReadProtected}
```

FullSphere is an option for GenerateGradients. If set True the gradients will be optimized on a full sphere rather than half a sphere. >>

```
Attributes[FullSphere] = {Protected, ReadProtected}
```

GradType is what type of gradient set will be produced in GenerateGradients "Normal" or "OverPlus". >>

```
Attributes[GradType] = {Protected, ReadProtected}
```

Method is an option for various algorithm-intensive functions that specifies what internal methods they should use. >>

```
Attributes[Method] = {Protected}
```

MethodReg is an options for RegisterData, RegisterDiffusionData, RegisterCardiacData and RegisterDataTransform. It specifies which registration method to use. Methods can be "rigid", "affine", "bspline" or "cyclic". >>

```
Attributes[MethodReg] = {Protected, ReadProtected}
```

OrderSpan is an options for FindOrder. >>

```
Attributes[OrderSpan] = {Protected, ReadProtected}
```

OutputPlot is an option for GradBmatrix. It specifies if the plots of the gradients should also be exported. >>

```
Attributes[OutputPlot] = {Protected, ReadProtected}
```

OutputType is an option for BmatrixCalc. Values can be "Matrix" or "Gradients". >>

```
Attributes[OutputType] = {Protected, ReadProtected}
```

PhaseEncoding is an options of GradSeq. Values can be "A", "P", "R" and "L". >>

```
Attributes[PhaseEncoding] = {Protected, ReadProtected}
```

Runs is an option for GenerateGradients. Set how often the minimalization function is run. The best solution of all runs is the output. Default value is 1. >>

```
Attributes[Runs] = {Protected, ReadProtected}
```

Steps is the number of step that is used in Generate Grads. >>

```
Attributes[Steps] = {Protected, ReadProtected}
```

StepSize is an option for GradBmatrix.

Specifies the integration stepsize is Method \rightarrow "Numerical" is used. \gg

```
Attributes[StepSize] = {Protected, ReadProtected}
```

SwitchAxes is an option for GradSeq. Default value is $\{\{1,2,3\},\{1,2,3\}\}$. First three

values are for diffusion gradients last three are for the acquisition gradients. \gg

```
Attributes[SwitchAxes] = {Protected, ReadProtected}
```

UnitMulti is an option for GradSeq. Default

value is 10^{-3} . Defines the scaling of the gradient strength. \gg

```
Attributes[UnitMulti] = {Protected, ReadProtected}
```

UseGrad is an option for GradSeq. The default value

is $\{0, 1, \{1, 0\}, 1\}$ where $\{\text{grex}, \text{gr180}, \{\text{grepi1}, \text{grepi2}\}, \text{grdiff}, \text{grflow}\}$. \gg

```
Attributes[UseGrad] = {Protected, ReadProtected}
```

VisualOpt is an option for GenerateGradients. Show

the minimalization proces of eacht calculation step. Default is False. \gg

```
Attributes[VisualOpt] = {Protected, ReadProtected}
```

ImportTools

Functions

BvalRead[file] imports the bvalue from a .dcm file. file must be a string. \gg

```
Attributes[BvalRead] = {Protected, ReadProtected}
```

```
SyntaxInformation[BvalRead] = {ArgumentsPattern  $\rightarrow$   $\{\_ \}$ }
```

GradRead[filename] imports the diffusion gradient direction from a .dcm file.

filename must be a string. \gg

```
Attributes[GradRead] = {Protected, ReadProtected}
```

```
Options[GradRead] = {ConvertDcm  $\rightarrow$  True}
```

```
SyntaxInformation[GradRead] = {ArgumentsPattern  $\rightarrow$   $\{\_, \text{OptionsPattern}[]\}$ }
```

ReadBrukerDiff[""] imports the bruker diffusion data selected by the input dialog.

ReadBrukerDiff["file"] imports the bruker diffusion data from "file", file must be location of 2dseq. \gg

```
Attributes[ReadBrukerDiff] = {Protected, ReadProtected}
```

```
Options[ReadBrukerDiff] = {BmatrixOut  $\rightarrow$  True}
```

```
SyntaxInformation[ReadBrukerDiff] = {ArgumentsPattern  $\rightarrow$   $\{\_, \text{OptionsPattern}[]\}$ }
```

`ReadBValue[folder,nr]` imports the gradient directions from the dicom header of the first nr of files in the given folder. folder must be a string, nr must be a int. Uses `BvalRead`. >>

```
Attributes[ReadBValue] = {Protected, ReadProtected}
```

```
SyntaxInformation[ReadBValue] = {ArgumentsPattern -> {_, _}}
```

`ReadDicom[folder]` imports all dicom files from the given folder.
`ReadDicom[{file1, file2,...}]` imports all the given filenames.
`ReadDicom[folder, {file1, file2,...}]` imports all the given filenames from the given folder.
`ReadDicom[folder, partsize]` imports all dicom files from the given folder and partitions them in given partsize.
`ReadDicom[{file1, file2, ...}, partsize]` imports all the given filenames and partitions them in given partsize.
`ReadDicom[folder, {file1, file2, ...}, partsize]` imports all the given filenames from the given folder and partitions them in given partsize. >>

```
Attributes[ReadDicom] = {Protected, ReadProtected}
```

```
Options[ReadDicom] = {ScaleCorrect -> False}
```

```
SyntaxInformation[ReadDicom] = {ArgumentsPattern -> {_, _., _., OptionsPattern[]}}
```

`ReadDicomDiff[folder, part]` imports all dicom files from the given folder and the corresponding diffusion parameters.

part is the number of diffusion images per slice including the unweighted images. >>

```
Attributes[ReadDicomDiff] = {Protected, ReadProtected}
```

```
Options[ReadDicomDiff] = {ScaleCorrect -> False}
```

```
SyntaxInformation[ReadDicomDiff] = {ArgumentsPattern -> {_, _, OptionsPattern[]}}
```

`ReadDicomDir[file]` reads the image data from a dicom directory. >>

```
Attributes[ReadDicomDir] = {Protected, ReadProtected}
```

```
SyntaxInformation[ReadDicomDir] = {ArgumentsPattern -> {_, _}}
```

`ReadDicomDirDiff[file]` reads the image data and relevant diffusion parameters from a dicom directory. >>

```
Attributes[ReadDicomDirDiff] = {Protected, ReadProtected}
```

```
Options[ReadDicomDirDiff] = {RotateGradient -> True}
```

```
SyntaxInformation[ReadDicomDirDiff] = {ArgumentsPattern -> {_, OptionsPattern[]}}
```


`ReadGradients[folder, nr]` imports the diffusion gradient directions from the dicom header of the first nr of files in the given folder.

folder must be a string, nr must be a int. Uses `GradRead`. >>

```
Attributes[ReadGradients] = {Protected, ReadProtected}
```

```
SyntaxInformation[ReadGradients] = {ArgumentsPattern -> {_, _}}
```

`ReadVoxSize[filename]` imports the voxelsize from a .dcm file. filename must be a string. Imports the pixel and slice spacing from the dicom header. Output is a list containing the voxels size {slice thickness, x, y}. >>

```
Attributes[ReadVoxSize] = {Protected, ReadProtected}
```

```
SyntaxInformation[ReadVoxSize] = {ArgumentsPattern -> {_}}
```

`ShiftPar[B0file.dcm, DTIfile.dcm]` imports the parameters from the dicom header and calculates the needed values to perform B0 field map correction. Needs a B0 dicom file and a diffusion dicom file. >>

```
Attributes[ShiftPar] = {Protected, ReadProtected}
```

```
SyntaxInformation[ShiftPar] = {ArgumentsPattern -> {_, _}}
```

Options

`BmatrixOut` is an option for `ImportBrukerData` if `True` the bmatrix is given, if `false` the gradients and `bvec` are given. >>

```
Attributes[BmatrixOut] = {Protected, ReadProtected}
```

`ConvertDcm` is an option for `GradRead`. >>

```
Attributes[ConvertDcm] = {Protected, ReadProtected}
```

`RotateGradient` is an option for `ReadDicomDirDiff`. If `False` it will also output the gradient direction as stored in the dicom header. >>

```
Attributes[RotateGradient] = {Protected, ReadProtected}
```

`ScaleCorrect` is an option for `ReadDicom`, `ReadDicomDiff`, `ReadDicomDir` and `ReadDicomDirDiff`. The dicom image values are corrected for rescale slope, scale slope and rescale intercept. >>

```
Attributes[ScaleCorrect] = {Protected, ReadProtected}
```

IVIMTools

Functions

BayesianIVIMFit2[data, bval, init, mask] performs bayesian IVIM fit of data.

data is the data which should be {slice, Ndiff, x, y}.

bval is the bvector whould be length Ndiff.

init is the initalization of the bayesian

fit which comes from IVIMCalc, (without S0 using 2 compartments).

mask is the region in which the bayesian fit is performed.

output is {f1, dc, pdc1}. >>

```
Attributes[BayesianIVIMFit2] = {Protected, ReadProtected}
```

```
Options[BayesianIVIMFit2] = {ChainSteps → {20000, 1000, 10}, UpdateStep → {0.5, 0.2, 0.5},
  FixPseudoDiff → False, CorrectPar → True, FixPseudoDiffSD → 0.5, OutputSamples → False,
  FitConstrains → ThetaConv[{{-7.6, 7.6}, {-10., -5.7}, {-7., 0.}}]}
```

```
SyntaxInformation[BayesianIVIMFit2] = {ArgumentsPattern → {_, _, _, _}, OptionsPattern[]}
```

BayesianIVIMFit3[data, bval, init, mask] performs bayesian IVIM fit of data.

data is the data which should be {slice, Ndiff, x, y}.

bval is the bvector whould be length Ndiff.

init is the initalization of the bayesian

fit which comes from IVIMCalc, (without S0 using 3 compartments).

mask is the region in which the bayesian fit is performed.

output is {f1, f2, dc, pdc1, pdc2}. >>

```
Attributes[BayesianIVIMFit3] = {Protected, ReadProtected}
```

```
Options[BayesianIVIMFit3] =
  {ChainSteps → {20000, 1000, 10}, UpdateStep → {0.5, 0.5, 0.1, 0.5, 0.5}, FixPseudoDiff → False,
  CorrectPar → True, OutputSamples → False, FixPseudoDiffSD → 0.5, FitConstrains →
  ThetaConv[{{-7.6, 7.6}, {-7.6, 7.6}, {-10., -5.5}, {-6.5, -2.3}, {-5.2, 0.}}]}
```

```
SyntaxInformation[BayesianIVIMFit3] = {ArgumentsPattern → {_, _, _, _}, OptionsPattern[]}
```

CorrectParMap[par, constraints, mask] removes

the IVIM parameters outside the constraints within the mask.

par is {f1, dc, pdc1} or {f1, f2, dc, pdc1, pdc2}.

constraints are the lower and upper constraints for each parameters {{min, max},...}

mask has the same dimensions as the parameter maps.

output are the corrected paremeter maps. >>

```
Attributes[CorrectParMap] = {Protected, ReadProtected}
```

```
SyntaxInformation[CorrectParMap] = {ArgumentsPattern → {_, _, _}}
```

FConvert[F] converts the fraction F from log space. >>

```
Attributes[FConvert] = {Protected, ReadProtected}
```

```
SyntaxInformation[FConvert] = {ArgumentsPattern → {_}}
```

FConverti[f] converts the fraction f to log space. >>

```
Attributes[FConverti] = {Protected, ReadProtected}
```

```
SyntaxInformation[FConverti] = {ArgumentsPattern → {_}}
```

FracCorrect[fraction, time] corrects the signal fraction calculated with the IVIM model for tissue relaxation and acquisition parameters.

After correction the signal fraction can be regarded as volume fraction.

FracCorrect[{fraction1, fraction2}, time] corrects the signal fraction1 and fraction2 from a 3 compartement IVIM model.

time is {{te, tr}, {t2t, t21}, {t1t, t11}} or {{te, tr}, {t2t, t21, t22}, {t1t, t11, t12}}

where t2t and t1t are "tissue"

relaxation times and t11 t12, t21 and t22 the "fluid" relaxation times

output is the corrected fraction maps >>

```
Attributes[FracCorrect] = {Protected, ReadProtected}
```

```
SyntaxInformation[FracCorrect] = {ArgumentsPattern → {_, _, _ .}}
```

HistogramPar[data, {constraints, Nbins}, style, color, range] plots histograms of IVIM solution.

HistogramPar[data, {constraints, Nbins, mu, conv}, components, color, range] plots histograms of IVIM solution.

data is {f1, dc, pdc1} or {f1, f2, dc, pdc1, pdc2}.

constraints are the ranges of the x-axes for the plots.

Nbins are the number of histogram bins.

style is the plot type, can be 1, 2, or 3.

color is the color of the histogram.

range are the ranges of the y-axes.

output is a row of histograms. >>

```
Attributes[HistogramPar] = {Protected, ReadProtected}
```

```
SyntaxInformation[HistogramPar] = {ArgumentsPattern → {_, _, _, _, _ .}}
```

IVIMCalc[data, binp, init] calculates the IVIM fit.

data should be 1D, 2D, 3D or 4D.

binp should be full bmatrix which can be calculated from the bvecs and bvals using Bmatrix.

init should be the initialization parameters for 2 components

this is {S0, f, D, Dp} for 2 components this is {S0, f1, f2, D, Dp1, Dp2}.

output is {S0, f1, D, pD1} or {S0, f1, f2, D, pD1, pD2}. >>

```
Attributes[IVIMCalc] = {Protected, ReadProtected}
```

```
Options[IVIMCalc] = {Method → Automatic, Parallelize → True, MonitorIVIMCalc → True,
  IVIMFixed → False, IVIMConstrained → True, IVIMTensFit → False, IVIMComponents → 2,
  IVIMConstrains → {{0.8, 1.2}, {0, 1}, {0.0005, 0.0035}, {0.001, 0.5}, {0.001, 0.5}}}
```

```
SyntaxInformation[IVIMCalc] = {ArgumentsPattern → {_, _, _, _}, OptionsPattern[]}
```

IVIMCorrectData[data, {S0, f, pdc}, bval] removes the ivim signal from the data.

data is the original data.

{S0, f, pdc} are the solution to a 2 compartment IVIM fit using IVIMCalc or BayesianIVIMFit2.

bval are the bvalues.

output is the corrected data. >>

```
Attributes[IVIMCorrectData] = {Protected, ReadProtected}
```

```
Options[IVIMCorrectData] = {FilterMaps → True, FilterType → Median, FilterSize → 1}
```

```
SyntaxInformation[IVIMCorrectData] = {ArgumentsPattern → {_, {_, _, _}, _, OptionsPattern[]}}
```

IVIMFunction[] gives the IVIM function with 2 comps.

IVIMFunction[components] gives the IVIM function.

IVIMFunction[components, type] gives the IVIM function.

type can be "Normal" or "Exp".

components can be 2 or 3.

output is the function with b, S0, f1, f2, D, pD1, pD2 as parameters >>

```
Attributes[IVIMFunction] = {Protected, ReadProtected}
```

```
SyntaxInformation[IVIMFunction] = {ArgumentsPattern → {_, _}}
```

IVIMResiduals[data, binp, pars] calculates the root mean square

residuals of an IVIM fit using IVIMCalc, BayesianIVIMFit2 or BayesianIVIMFit3. >>

```
Attributes[IVIMResiduals] = {Protected, ReadProtected}
```

ThetaConv[{F1, Fc, pDc}] converts the parameters from Log space to normal space. Is used in BayesianIVIMFit2 and BayesianIVIMFit3.

ThetaConv[{F1, F2, Dc, pDc1}] converts the parameters from Log space to normal space. Is used in BayesianIVIMFit2 and BayesianIVIMFit3.

ThetaConv[{F1, F2, Dc, pDc1, pDc2}] converts the parameters from Log space to normal space. Is used in BayesianIVIMFit2 and BayesianIVIMFit3. >>

Attributes[ThetaConv] = {Protected, ReadProtected}

ThetaConvi[{f, dc, pdc}] converts the parameters from Normal space to Log space. Is used in BayesianIVIMFit2 and BayesianIVIMFit3.

ThetaConvi[{f1, f2, dc, pdc1}] converts the parameters from Normal space to Log space. Is used in BayesianIVIMFit2 and BayesianIVIMFit3.

ThetaConvi[{f1, f2, dc, pdc1, pdc2}] converts the parameters from Normal space to Log space. Is used in BayesianIVIMFit2 and BayesianIVIMFit3. >>

Attributes[ThetaConvi] = {Protected, ReadProtected}

SyntaxInformation[ThetaConvi] = {ArgumentsPattern → {_}}

Options

ChainSteps is an option for BayesianIVIMFit2 and BayesianIVIMFit3. It determines how long the algorithm runs. three values must be given {iterations, burn steps, sample density}. >>

Attributes[ChainSteps] = {Protected, ReadProtected}

CorrectPar is an option for BayesianIVIMFit2 and BayesianIVIMFit3. If True it removes the values outside the constraints using CorrectParMap >>

Attributes[CorrectPar] = {Protected, ReadProtected}

FilterMaps is an option for IVIMCorrectData. If True the IVIM parameter maps are filtered before signal correction >>

Attributes[FilterMaps] = {Protected, ReadProtected}

FilterSize is an option for IVIMCorrectData. If FilterMaps is True it gives the kernel size. >>

Attributes[FilterSize] = {Protected, ReadProtected}

FilterType is an option for IVIMCorrectData. If FilterMaps is True it tells which filter to use. can be "Median" of "Gaussian" >>

Attributes[FilterType] = {Protected, ReadProtected}

FitConstrains is an option for BayesianIVIMFit2 and BayesianIVIMFit3. Gives the constraints of the parameters. The values are used for displaying the histograms and for the initialization if CorrectPar is True >>

Attributes[FitConstrains] = {Protected, ReadProtected}

FixPseudoDiff is an option for BayesianIVIMFit2 and BayesianIVIMFit3. If the pDc1 and pD2 were fixed in IVIMCalc this value should be True. >>

Attributes[FixPseudoDiff] = {Protected, ReadProtected}

FixPseudoDiffSD is an option for BayesianIVIMFit2 and BayesianIVIMFit3. Gives the standard deviation of pDc1 and pD2 if FixPseudoDiff is True >>

Attributes[FixPseudoDiffSD] = {Protected, ReadProtected}

IVIMComponents is an option for IVIMCalc. Default value is 2, the tissue and the blood component. can also be set to 3. >>

Attributes[IVIMComponents] = {Protected, ReadProtected}

IVIMConstrained is an option for IVIMCalc. When set True the fit will be constrained to the values given in IVIMConstrains. >>

Attributes[IVIMConstrained] = {Protected, ReadProtected}

IVIMConstrains is an option for IVIMCalc. Default values are: {{0.8, 1.2}, {0, 1}, {0.0005, 0.0035}, {0.005, 0.5}, {0.002, 0.015}}. Where {{S0 in percentage},{fractions},{tissue diffusion},{blood compartment Dp},{third compartment}}. >>

Attributes[IVIMConstrains] = {Protected, ReadProtected}

IVIMFixed is an option for IVIMCalc and the default value is False. When set True the pseudo diffusion will be fixed to the parameter given as init. When set to "One" only the fast component of a 3 compartment fit is fixed. >>

Attributes[IVIMFixed] = {Protected, ReadProtected}

IVIMTensFit is an option for IVIMCalc. When set True the tissue diffusion component will be calculated as a tensor. >>

Attributes[IVIMTensFit] = {Protected, ReadProtected}

Method is an option for various algorithm-intensive functions that specifies what internal methods they should use. >>

Attributes[Method] = {Protected}

MonitorIVIMCalc is an option for IVIMCalc. When true the proceses of the calculation is shown. >>

```
Attributes[MonitorIVIMCalc] = {Protected, ReadProtected}
```

OutputSamples is an option for BayesianIVIMFit2 and BayesianIVIMFit3.

If set True the full marcov chain is given as an additional output. >>

```
Attributes[OutputSamples] = {Protected, ReadProtected}
```

Parallelize[*expr*] evaluates *expr* using automatic parallelization. >>

```
Attributes[Parallelize] = {HoldFirst, Protected}
```

```
Parallelize[Parallel`Kernels`Private`args$___] :=  
  (Parallel`Protected`doAutoLaunch[TrueQ[Parallel`Static`$enableLaunchFeedback]];  
   Parallelize[Parallel`Kernels`Private`args$])
```

```
Options[Parallelize] = {Method → Automatic, DistributedContexts → $Context}
```

UpdateStep is an option for BayesianIVIMFit2 and BayesianIVIMFit3. It determines how often the parameters are updated. Is optimized during the first 500 burn steps. >>

```
Attributes[UpdateStep] = {Protected, ReadProtected}
```

JcouplingTools

Functions

GetSpinSystem[name] get a spinsystem that can be used in SimHamiltonian. Current implementes systems are "glu", "lac", "gaba", "fatGly", "fatAll", "fatEnd", "fatDouble", "fatStart", and "fatMet". >>

```
Attributes[GetSpinSystem] = {Protected, ReadProtected}
```

```
Options[GetSpinSystem] = {CenterFrequency → 4.65}
```

```
SyntaxInformation[GetSpinSystem] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

PhaseAlign[spec] automatically phase aligns the spectrum by maximizing the Real part of the spectrum. >>

```
Attributes[PhaseAlign] = {Protected, ReadProtected}
```

```
SyntaxInformation[PhaseAlign] = {ArgumentsPattern → {_}}
```

PlotSpectrum[ppm, spec] plots the spectrum, ppm and spec can be generated using SimReadout. >>

```
Attributes[PlotSpectrum] = {Protected, ReadProtected}
```

```
Options[PlotSpectrum] = {PlotRange → {{0, 6}, Full}, SpectrumColor → GrayLevel[0]}
```

```
SyntaxInformation[PlotSpectrum] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

SequencePulseAcquire[din, H] performs a pulsaquire experiment of the spin system din given the hamiltonian H with a 90 Degree pulse.
 SequencePulseAcquire[din, H, b1] performs a pulsaquire experiment of the spin system din given the hamiltonian H with a 90 Degree pulse and b1.
 The output is a new spinsystem dout. >>

```
Attributes[SequencePulseAcquire] = {Protected, ReadProtected}
```

```
SyntaxInformation[SequencePulseAcquire] = {ArgumentsPattern -> {_, _, _..}}
```

SequenceSpinEcho[din, H, te] performs a spin echo experiment with echo time te of the spin system din given the hamiltonian H with a 90 and 180 Degree pulse.
 SequenceSpinEcho[din, H, te, b1] performs a spin echo experiment with echo time te of the spin system din given the hamiltonian H with a 90 and 180 Degree pulse and b1.
 The output is a new spinsystem dout. >>

```
Attributes[SequenceSpinEcho] = {Protected, ReadProtected}
```

```
SyntaxInformation[SequenceSpinEcho] = {ArgumentsPattern -> {_, _, _, _..}}
```

SequenceSteam[din, H, {te, tm}] performs a stimulated echo experiment with echo time te and mixing time tm of the spin system din given the hamiltonian H with 3 90 Degree pulses.
 The output is a new spinsystem dout. >>

```
Attributes[SequenceSteam] = {Protected, ReadProtected}
```

```
SyntaxInformation[SequenceSteam] = {ArgumentsPattern -> {_, _, {_, _}}}
```

SequenceTSE[din, H, {te, necho}, {ex, ref}] performs a multi echo spin echo experiment with echo time te with necho echos of the spin system din given the hamiltonian H using ex Degree exitation and ref Degree refocus pulses.
 SequenceTSE[din, H, {te, necho}, {ex, ref}, b1_:1] performs a multi echo spin echo experiment with echo time te with necho echos of the spin system din given the hamiltonian H using ex Degree exitation and ref Degree refocus pulses and b1.
 The output is a new spinsystem dout. >>

```
Attributes[SequenceTSE] = {Protected, ReadProtected}
```

```
SyntaxInformation[SequenceTSE] = {ArgumentsPattern -> {_, _, {_, _}, {_, _}, _..}}
```

SimAddPhase[din, H, phase] adds phase to the spin system din given the hamiltonian H. din and H are generated by SimHamiltonian.
 The output is a new spinsystem dout. >>

```
Attributes[SimAddPhase] = {Protected, ReadProtected}
```

```
SyntaxInformation[SimAddPhase] = {ArgumentsPattern -> {_, _, _}}
```


`SimEvolve[din,H,t]` evolves the spin system `din` given the hamiltonian `H` over a time `t`. `din` and `H` are generated by `SimHamiltonian`. The output is a new spinsystem `dout`. >>

```
Attributes[SimEvolve] = {Protected, ReadProtected}
```

```
SyntaxInformation[SimEvolve] = {ArgumentsPattern -> {_, _, _}}
```

`SimHamiltonian[sysi]` simulates the hamiltonian for a given spin system. The spinsystem is generated by `GetSpinSystem`. The output is the spin system and hamiltonian structure. >>

```
Attributes[SimHamiltonian] = {Protected, ReadProtected}
```

```
Options[SimHamiltonian] = {FieldStrength -> 3}
```

```
SyntaxInformation[SimHamiltonian] = {ArgumentsPattern -> {_, OptionsPattern[]}}
```

`SimReadout[din, H]` performs a readout of a spinsystem `din` with hamiltonian `H`. Output is `{time,fids,ppm,spec,dout}`, which are the free induction decay `fids` with its time, the spectrum `spec` with its ppm and the evolved spin system `dout`. >>

```
Attributes[SimReadout] = {Protected, ReadProtected}
```

```
Options[SimReadout] = {ReadoutOutput -> all, ReadoutPhase -> 90,  
  Linewidth -> 5, LinewidthShape -> L, ReadoutSamples -> 2046, ReadoutBandwidth -> 2000}
```

```
SyntaxInformation[SimReadout] = {ArgumentsPattern -> {_, _, OptionsPattern[]}}
```

`SimRotate[din, H ,angle]` rotates the spin system `din` given the hamiltonian `H` over `angle` with phase 90 degrees. `SimRotate[din, H ,angle, phase]` rotates the spin system `din` given the hamiltonian `H` over `angle` with `phase`. `din` and `H` are generated by `SimHamiltonian`. The output is a new spinsystem `dout`. >>

```
Attributes[SimRotate] = {Protected, ReadProtected}
```

```
SyntaxInformation[SimRotate] = {ArgumentsPattern -> {_, _, _, _}}
```

`SimSignal[din, H]` performs a readout of a spinsystem `din` with hamiltonian `H`. Output is the complex signal. >>

```
Attributes[SimSignal] = {Protected, ReadProtected}
```

```
Options[SimSignal] = {ReadoutOutput -> all}
```

```
SyntaxInformation[SimSignal] = {ArgumentsPattern -> {_, _, OptionsPattern[]}}
```

SimSpoil[din] spoils all the non zeroth order states of a spin system.
The output is a new spinsystem dout. >>

```
Attributes[SimSpoil] = {Protected, ReadProtected}
```

```
SyntaxInformation[SimSpoil] = {ArgumentsPattern -> { _ } }
```

SysTable[sys] shows the spinsystem as a table. The spinsytem is obtained form GetSpinSystem. >>

```
Attributes[SysTable] = {Protected, ReadProtected}
```

```
SyntaxInformation[SysTable] = {ArgumentsPattern -> { _ } }
```

Options

CenterFrequency is an option for GetSpinSystem and defines the center frequency. >>

```
Attributes[CenterFrequency] = {Protected, ReadProtected}
```

FieldStrength is an option for SimHamiltonian. It
defines the field strength for which the hamiltonian is calculated. >>

```
Attributes[FieldStrength] = {Protected, ReadProtected}
```

Linewidth is an option for SimReadout and defines the spectral linewidth. >>

```
Attributes[Linewidth] = {Protected, ReadProtected}
```

LinewidthShape is an option for SimReadout and defines the linewidth shape, values can be
"L", "G" or "L", which are Laplacian, Gaussian or a combination, respectively. >>

```
Attributes[LinewidthShape] = {Protected, ReadProtected}
```

PlotRange is an option for graphics functions
that specifies what range of coordinates to include in a plot. >>

```
Attributes[PlotRange] = {Protected, ReadProtected}
```

ReadoutBandwith is an option for SimReadout defines the spectral bandwith. >>

```
Attributes[ReadoutBandwith] = {Protected, ReadProtected}
```

ReadoutOutput is an option for SimReadout and SimSignal and values
can be "all" and "each". When set to "all" the total signal and signal is given,
when set to "each" the signal or spectrum for each peak is given seperately. >>

```
Attributes[ReadoutOutput] = {Protected, ReadProtected}
```

ReadoutPhase is an option for SimReadout and defines the readout phase. >>

```
Attributes[ReadoutPhase] = {Protected, ReadProtected}
```

`ReadoutSamples` is an option for `SimReadout` and defines the number of readout samples for the spectrum. >>

```
Attributes[ReadoutSamples] = {Protected, ReadProtected}
```

`SpectrumColor` is an option for `PlotSpectrum` and defines the spectrum color. >>

```
Attributes[SpectrumColor] = {Protected, ReadProtected}
```

MaskingTools

Functions

`GetMaskData[data, mask]` retruns the data selected by the mask. >>

```
Attributes[GetMaskData] = {Protected, ReadProtected}
```

```
Options[GetMaskData] = {GetMaskOutput → All}
```

```
SyntaxInformation[GetMaskData] = {ArgumentsPattern → {_, _}, OptionsPattern[]}
```

`HomoginizeData[data, mask]` tries to homoginize the data within the mask by removing intensity gradients. >>

```
Attributes[HomoginizeData] = {Protected, ReadProtected}
```

```
SyntaxInformation[HomoginizeData] = {ArgumentsPattern → {_, _}}
```

`Mask[data]` creates a mask by automatically finding a threshold.

`Mask[data, min]` creates a mask which selects only data above the min value.

`Mask[data, {min, max}]` creates a mask which selects data between the min and max value. >>

```
Attributes[Mask] = {Protected, ReadProtected}
```

```
Options[Mask] = {MaskSmoothing → False, MaskComponents → 1, MaskClosing → 5, MaskFiltKernel → 2}
```

```
SyntaxInformation[Mask] = {ArgumentsPattern → {_, _}, OptionsPattern[]}
```

`MaskData[data, mask]` applies a mask to data. mask can be 2D or 3D, data can be 2D, 3D or 4D. >>

```
Attributes[MaskData] = {Protected, ReadProtected}
```

```
SyntaxInformation[MaskData] = {ArgumentsPattern → {_, _}}
```

`MeanSignal[data]` calculates the mean signal per volume of 4D data.

`MeanSignal[data, pos]` calculates the

mean signal per volume of 4D data for the given list of positions. >>

```
Attributes[MeanSignal] = {Protected, ReadProtected}
```

```
Options[MeanSignal] = {UseMask → True}
```

```
SyntaxInformation[MeanSignal] = {ArgumentsPattern → {_, _}, OptionsPattern[]}
```

`MergeSegmentations[masks, labels]` generates an ITKsnap or slices3D compatible segmentation from individual masks and label numbers. Output is a labeled segmentation. >>

`Attributes[MergeSegmentations] = {Protected, ReadProtected}`

`SyntaxInformation[MergeSegmentations] = {ArgumentsPattern → {_, _}}`

`NormalizeData[data]` normalizes the data to the mean signal of the data. For 4D data it normalizes to the first volume of the 4th dimension. `NormalizeData[data, {min, max}]` normalizes the data between min and max. >>

`Attributes[NormalizeData] = {Protected, ReadProtected}`

`SyntaxInformation[NormalizeData] = {ArgumentsPattern → {_, _., OptionsPattern[]}}`

`RemoveMaskOverlaps[mask]` removes the overlaps between multiple masks. Mask is a 4D dataset with {z, masks, x, y} >>

`Attributes[RemoveMaskOverlaps] = {Protected, ReadProtected}`

`SyntaxInformation[RemoveMaskOverlaps] = {ArgumentsPattern → {_}}`

`RescaleSegmentation[data, dim]` rescales segmentations to given dimensions. `RescaleSegmentation[data, {vox1, vox2}]` rescales segmentations from voxelsize vox1 to voxelsize vox2. >>

`Attributes[RescaleSegmentation] = {Protected, ReadProtected}`

`SyntaxInformation[RescaleSegmentation] = {ArgumentsPattern → {_, _}}`

`ROIMask[maskdim, {name→{{x,y},slice}..}]` crates mask from coordinates x and y at slice. maskdim is the dimensions of the output {zout,xout,yout}. >>

`Attributes[ROIMask] = {Protected, ReadProtected}`

`SyntaxInformation[ROIMask] = {ArgumentsPattern → {_, _, _}}`

`SegmentMask[mask, n]` devides a mask in n equal segments along the slice direction. n must be an integer. >>

`Attributes[SegmentMask] = {Protected, ReadProtected}`

`SyntaxInformation[SegmentMask] = {ArgumentsPattern → {_, _, _}}`

`SmoothMask[mask]` generates one clean masked volume form a noisy mask. >>

`Attributes[SmoothMask] = {Protected, ReadProtected}`

`Options[SmoothMask] = {MaskComponents → 1, MaskClosing → 5, MaskFiltKernel → 2}`

`SyntaxInformation[SmoothMask] = {ArgumentsPattern → {_, OptionsPattern[]}}`

SmoothSegmentation[masks] smooths segmentations
and removes the overlaps between multiple segmentations. >>

Attributes[SmoothSegmentation] = {Protected, ReadProtected}

Options[SmoothSegmentation] = {MaskFiltKernel → 2}

SyntaxInformation[SmoothSegmentation] = {ArgumentsPattern → {_, OptionsPattern[]}}

SplitSegmentations[segmentation] splits a lable
mask from ITKsnap or slicer3D in separte masks and label numbers.
Output is masks and label numbers, {mask, labs}. >>

Attributes[SplitSegmentations] = {Protected, ReadProtected}

SyntaxInformation[SplitSegmentations] = {ArgumentsPattern → {_}}

Options

GetMaskOutput is an option for GetMaskData. Defaul is "Slices" which
gives the mask data per slices. Else the entire mask data is given as output. >>

Attributes[GetMaskOutput] = {Protected, ReadProtected}

MaskClosing is an option for Mask and
SmoothMask. The size of the holes in the mask that will be closed >>

Attributes[MaskClosing] = {Protected, ReadProtected}

MaskComponents is an option for Mask and
SmoothMask. Determinse the amount of largest clusters used as mask. >>

Attributes[MaskComponents] = {Protected, ReadProtected}

MaskFiltKernel is an option for Mask, SmoothMask
and SmoothSegmentation. How mucht the contours are smoothed. >>

Attributes[MaskFiltKernel] = {Protected, ReadProtected}

MaskSmoothing is an options for Mask, if set to True
it smooths the mask, by closing holse and smoothing the contours. >>

Attributes[MaskSmoothing] = {Protected, ReadProtected}

UseMask is a function for MeanSignal and DriftCorrect >>

Attributes[UseMask] = {Protected, ReadProtected}

NiftiTools

Functions

CompressNiiFiles[] prompts for a folder. It
 then compresses all nii files to .nii.gz files in the selected folder.
 CompressNiiFiles[folder] compresses all nii files to .nii.gz files in folder. >>

```
Attributes[CompressNiiFiles] = {Protected, ReadProtected}
```

```
SyntaxInformation[CompressNiiFiles] = {ArgumentsPattern → {_, _}}
```

DcmToNii[] converts a dicom folder to nii, you will be promoted for the location of the folders.
 DcmToNii[{"input", "output"}] converts the
 "input" dicom folder to nii files which are place in the "output" folder.
 For this function to work the dcm2niix.exe file should be present in the QMRITools application folder. >>

```
Attributes[DcmToNii] = {Protected, ReadProtected}
```

```
Options[DcmToNii] = {CompressNii → True, Method → Automatic}
```

```
SyntaxInformation[DcmToNii] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

ExportBmat[bmat] exports the diffusion bmatrix to exploreDTI format.
 ExportBmat[bmat, "file"] exports the diffusion bmatrix to "file" in the exploreDTI format. >>

```
Attributes[ExportBmat] = {Protected, ReadProtected}
```

```
SyntaxInformation[ExportBmat] = {ArgumentsPattern → {_, _}}
```

ExportBval[bvals] exports the diffusion bvalues to exploreDTI format.
 ExportBval[bvals, "file"] exports the diffusion bvalues to "file" in the exploreDTI format. >>

```
Attributes[ExportBval] = {Protected, ReadProtected}
```

```
SyntaxInformation[ExportBval] = {ArgumentsPattern → {_, _}}
```

ExportBvec[grad] exports the diffusion gradients to exploreDTI format.
 ExportBvec[grad, "file"] exports the diffusion gradients to "file" in the exploreDTI format. >>

```
Attributes[ExportBvec] = {Protected, ReadProtected}
```

```
SyntaxInformation[ExportBvec] = {ArgumentsPattern → {_, _}}
```

ExportNii[data, vox] exports the nii file and will prompt for a file name.
 ExportNii[data, vox, "file"] exports the nii file to the location "file". >>

```
Attributes[ExportNii] = {Protected, ReadProtected}
```

```
Options[ExportNii] = {NiiDataType → Automatic, CompressNii → True}
```

```
SyntaxInformation[ExportNii] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

ExtractNiiFiles[] prompts for a folder. It then extracts all nii.gz files to .nii files in the selected folder.
 ExtractNiiFiles[folder] extracts all nii.gz files to .nii files in folder. >>

```
Attributes[ExtractNiiFiles] = {Protected, ReadProtected}
```

```
SyntaxInformation[ExtractNiiFiles] = {ArgumentsPattern → {_., _.
```

ImportBmat[] will prompt to select the *.txt file containing the bmatrix.
 ImportBmat[*.txt] imports the given *.txt file containing the bmatrix. >>

```
Attributes[ImportBmat] = {Protected, ReadProtected}
```

```
SyntaxInformation[ImportBmat] = {ArgumentsPattern → {_., _.
```

ImportBval[] will prompt to select the *.bval file.
 ImportBval[*.bval] imports the given *.bval file. >>

```
Attributes[ImportBval] = {Protected, ReadProtected}
```

```
SyntaxInformation[ImportBval] = {ArgumentsPattern → {_., _.
```

ImportBvalvec[] will prompt to select the *.bval and *.bvec files.
 ImportBvalvec[file] if file is either a *.bval
 or *.bvec it will automatically import the *.bval and *.bvec files.
 ImportBvalvec[*.bvec,*.bval] imports the given *.bval and *.bvec files. >>

```
Attributes[ImportBvalvec] = {Protected, ReadProtected}
```

```
Options[ImportBvalvec] = {FlipBvec → False}
```

```
SyntaxInformation[ImportBvalvec] = {ArgumentsPattern → {_., _., OptionsPattern[]}}
```

ImportBvec[] will prompt to select the *.bvec file.
 ImportBvec[*.bvec] imports the given *.bvec file. >>

```
Attributes[ImportBvec] = {Protected, ReadProtected}
```

```
Options[ImportBvec] = {FlipBvec → False}
```

```
SyntaxInformation[ImportBvec] = {ArgumentsPattern → {_., OptionsPattern[]}}
```

ImportExploreDTIens["file"] imports the *.nii export for the tensor from explore DTI. >>

```
Attributes[ImportExploreDTIens] = {Protected, ReadProtected}
```

ImportNii[] prompts to select the nii file to import.
 ImportNii["file"] imports the nii file.
 The default output is {data, vox}, however using NiiMethod various outputs can be given.
 The Nii import is also supported using the native Import function from Mathematica. >>

```
Attributes[ImportNii] = {Protected, ReadProtected}
```

```
Options[ImportNii] = {NiiMethod → default, NiiScaling → False}
```

```
SyntaxInformation[ImportNii] = {ArgumentsPattern → {_., OptionsPattern[]}}
```

ImportNiiDiff[] will prompt for the *.nii, *.bvec and *.bval file to import.
 ImportNiiDiff[*.nii] will import the *.nii file and
 automatically also imports the *.bvec and *.bval if they have the same name.
 ImportNiiDiff[*.nii,*.bvec,*.bval] will import the given files.
 The output will be {data,grad,bvec,vox}. >>

```
Attributes[ImportNiiDiff] = {Protected, ReadProtected}
```

```
Options[ImportNiiDiff] = {RotateGradients → False, FlipBvec → True}
```

```
SyntaxInformation[ImportNiiDiff] = {ArgumentsPattern → {_., _., _., OptionsPattern[]}}
```

ImportNiiDix["file"] imports the dixon nii file which should contain
 all possible outputs given by the scanner and corrects them accordingly. >>

```
Attributes[ImportNiiDix] = {Protected, ReadProtected}
```

ImportNiiT1["file"] imports the T1 file which should contain the echos
 and the T1map calculated by the scanner and corrects them accordingly. >>

```
Attributes[ImportNiiT1] = {Protected, ReadProtected}
```

ImportNiiT2["file"] imports the T2 file which should contain the echos
 and the T2map calculated by the scanner and corrects them accordingly. >>

```
Attributes[ImportNiiT2] = {Protected, ReadProtected}
```

Options

CompressNii is an option for DcmToNii and ExportNii. If set True .nii.gz files will be created. >>

```
Attributes[CompressNii] = {Protected, ReadProtected}
```

FlipBvec is an option for ImportBvalvec. >>

```
Attributes[FlipBvec] = {Protected, ReadProtected}
```

Method is an option for various algorithm-intensive
 functions that specifies what internal methods they should use. >>

```
Attributes[Method] = {Protected}
```

NiiDataType is an option of Export Nii. The number
 type of Nii file can be "Integer", "Real", "Complex", or "Automatic". >>

```
Attributes[NiiDataType] = {Protected, ReadProtected}
```

NiiMethod is an option for ImportNii. Values can be
 "data", "dataTR", "header", "scaling", "headerMat", "rotation", "all". >>

```
Attributes[NiiMethod] = {Protected, ReadProtected}
```


NiiScaling is an option for ImportNii. It scales the nii values with scale slope and offset for quantitative data. >>

```
Attributes[NiiScaling] = {Protected, ReadProtected}
```

RotateGradients is an option for ImportNiiDiff. >>

```
Attributes[RotateGradients] = {Protected, ReadProtected}
```

PhysiologyTools

Functions

AlignRespLog[physLog, respirect, scanTime] aligns respirect and physlog data. physLog is output from ImportPhyslog. respirect is the first output from ImportRespirect. >>

```
Attributes[AlignRespLog] = {Protected, ReadProtected}
```

```
Options[AlignRespLog] = {OutputMethod → val, SampleStep → 0.005}
```

```
SyntaxInformation[AlignRespLog] = {ArgumentsPattern → {_, _, _}, OptionsPattern[]}
```

ImportPhyslog[] imports all physlog files from the folder selcted.

ImportPhyslog["folder"] imports all physlog files from "folder" selcted. >>

```
Attributes[ImportPhyslog] = {Protected, ReadProtected}
```

```
SyntaxInformation[ImportPhyslog] = {ArgumentsPattern → {_.}}
```

ImportRespirect[] impors all the respirect log files from the folder selcted.

ImportRespirect["folder"] impors all the respirect log files from the "folder" selcted. >>

```
Attributes[ImportRespirect] = {Protected, ReadProtected}
```

```
SyntaxInformation[ImportRespirect] = {ArgumentsPattern → {_.}}
```

PlotPhyslog[{time, resp}, {start, stop}] plots the physlog from ImportPhyslog.

PlotPhyslog[{time, resp}, {start, stop}, scanTime] plots the physlog from ImportPhyslog. >>

```
Attributes[PlotPhyslog] = {Protected, ReadProtected}
```

```
SyntaxInformation[PlotPhyslog] = {ArgumentsPattern → {_, _, _}}
```

PlotRespiract[data, dataP, scantimes] plots the respirect data to correct peaks. data and dataP are the first outputs of ImportResirect. scantimes is the output from AlignRespLog.

PlotRespiract[data, dataP, scantimes, steps]. >>

```
Attributes[PlotRespiract] = {Protected, ReadProtected}
```

```
SyntaxInformation[PlotRespiract] = {ArgumentsPattern → {_, _, _}}
```

Options

OutputMethod can be "val" or "plot" >>

```
Attributes[OutputMethod] = {Protected, ReadProtected}
```

SampleStep is an option for AlignRespiract >>

```
Attributes[SampleStep] = {Protected, ReadProtected}
```

PlottingTools

Functions

GetSliceData[data, offsets] gets the slices from
the data defined by offsets which are obtained by GetSlicePosisions.
GetSliceData[data, offsets, vox] gets the slices from the data defined by
offsets which are obtained by GetSlicePosisions in mm.
The offsets can also be provided manually which is {{AX,..},{COR,..},{SAG,..}}. >>

```
Attributes[GetSliceData] = {Protected, ReadProtected}
```

```
SyntaxInformation[GetSliceData] = {ArgumentsPattern -> {_, _, _..}}
```

GetSlicePositions[data] finds the position of slices with the maximal signal in voxel index.
GetSlicePositions[data, vox] find the position of slices with the maximal signal in mm. >>

```
Attributes[GetSlicePositions] = {Protected, ReadProtected}
```

```
Options[GetSlicePositions] =  
{MakeCheckPlot -> False, DropSlices -> {1, 1, 1}, PeakNumber -> {1, 1, 2}}
```

```
SyntaxInformation[GetSlicePositions] = {ArgumentsPattern -> {_, .., OptionsPattern[]}}
```

GradientPlot[bvec, bval] plots the given
bvec with position of the gradients scaled according to the bval. >>

```
Attributes[GradientPlot] = {Protected, ReadProtected}
```

```
Options[GradientPlot] =  
{PlotSpace -> bspace, PlotColor -> Auto, SphereSize -> 0.05, PositiveZ -> False}
```

```
SyntaxInformation[GradientPlot] = {ArgumentsPattern -> {_, .., OptionsPattern[]}}
```

ListSpherePlot[points] plots 3D points as spheres >>

```
Attributes[ListSpherePlot] = {Protected, ReadProtected}
```

```
Options[ListSpherePlot] = {SphereSize -> 2, SphereColor -> Automatic}
```

```
SyntaxInformation[ListSpherePlot] = {ArgumentsPattern -> {_, OptionsPattern[]}}
```

MakeSliceImages[imgData] generates images from the imgData which is obtained from GetSliceData.

MakeSliceImages[imgData, vox] generates images from the imgData which is obtained from GetSliceData, vox is used for the correct aspect ratio of the images.

MakeSliceImages[imgData, {labData, labels}] generates images from the imgData which is obtained from GetSliceData with an overlay of the segmentations in labData, which can also be obtained using GetSliceData on the segmentations.

labels should be the label numbers used in the original segmentation (to allow correct scaling between slices).

MakeSliceImages[imgData, {labData, labels}, vox] generates images from the imgData which is obtained from GetSliceData with an overlay of the segmentations in labData, which can also be obtained using GetSliceData on the segmentations, vox is used for the correct aspect ratio of the images. >>

```
Attributes[MakeSliceImages] = {Protected, ReadProtected}
```

```
Options[MakeSliceImages] =  
{PlotRange → Automatic, ColorFunction → GrayTones, ImageLegend → False}
```

```
SyntaxInformation[MakeSliceImages] = {ArgumentsPattern → {_, __, __, OptionsPattern[]}}
```

PlotContour[data, vox] creates a contour of the data.

PlotContour[data, vox, scale] creates a contour of the data with the surface colored according to scale.

PlotContour[data, vox, scale, range] creates a contour of the data with the surface colored according to scale with a fixed plotrange. >>

```
Attributes[PlotContour] = {Protected, ReadProtected}
```

```
Options[PlotContour] = {ContourStyle → {GrayLevel[0.5], 0.25}}
```

```
SyntaxInformation[PlotContour] = {ArgumentsPattern → {_, __, __, __, OptionsPattern[]}}
```

PlotCorrection[w] plots deformation vectors w {w1,w2..} generated by Registration2D and Registration3D for multiple datasets or registration steps. >>

```
Attributes[PlotCorrection] = {Protected, ReadProtected}
```

```
SyntaxInformation[PlotCorrection] = {ArgumentsPattern → {_}}
```

PlotData[data] plots the data.

PlotData[data, vox] plots the data and for 3D and 4D data assumes the voxelsize vox (z,x,y).

PlotData[data1, data2] plots data1 and data2.

PlotData[data1, data2, vox] plots data1 and data2 and for 3D and 4D data assumes the voxelsize vox (z,x,y). >>

```
Attributes[PlotData] = {Protected, ReadProtected}
```

```
Options[PlotData] = {PlotRange → Auto, ColorFunction → BlackToWhite}
```

```
SyntaxInformation[PlotData] = {ArgumentsPattern → {_, _., _., OptionsPattern[]}}
```

PlotData3D[data,vox] is a 3D dataviewer, data
is the 3D dataset and voxsize the size of the voxels in mm (z,x,y). >>

```
Attributes[PlotData3D] = {Protected, ReadProtected}
```

```
SyntaxInformation[PlotData3D] = {ArgumentsPattern → {_, _., _}}
```

PlotDefGrid[data, phasemap, shiftpar] plots the dataset on the background
with on top the non deformed and the deformed grid, or arrows or lines. >>

```
Attributes[PlotDefGrid] = {Protected, ReadProtected}
```

PlotDuty[{grad, bval, ord}, mode] plot the gradient dutycycle >>

```
Attributes[PlotDuty] = {Protected, ReadProtected}
```

```
SyntaxInformation[PlotDuty] = {ArgumentsPattern → {{_, _, _}, _., _}}
```

PlotIVIM[vals, data, bvals] plots the results of the IVIM fits from IVIMCalc or BayesianIVIMFit2 or Baye. >>

```
Attributes[PlotIVIM] = {Protected, ReadProtected}
```

```
Options[PlotIVIM] =  
{Method → , PlotColor → {RGBColor[1, 0, 0], RGBColor[0, 1, 0], RGBColor[0, 0, 1], GrayLevel[0]},  
NormalizeIVIM → Fit, PlotRange → Auto, ImageSize → 400}
```

```
SyntaxInformation[PlotIVIM] = {ArgumentsPattern → {_, _, _, OptionsPattern[]}}
```

PlotMoments[{G(t),...}, te, t] plots the moments generated by CalculateMoments >>

```
Attributes[PlotMoments] = {Protected, ReadProtected}
```

```
SyntaxInformation[PlotMoments] = {ArgumentsPattern → {_, _, _}}
```

PlotSequence[seq,var] where seq is the output from GradSeq. >>

```
Attributes[PlotSequence] = {Protected, ReadProtected}
```

```
SyntaxInformation[PlotSequence] = {ArgumentsPattern → {_, _}}
```

Options

ColorFunction is an option for graphics functions
that specifies a function to apply to determine colors of elements. >>

```
Attributes[ColorFunction] = {Protected}
```

ContourStyle is an option for contour plots that
specifies the style in which contour lines or surfaces should be drawn. >>

```
Attributes[ContourStyle] = {Protected}
```

DropSlices is an option for GetSlicePositions and specifies how many slices from the beginning and should be ignored. >>

```
Attributes[DropSlices] = {Protected, ReadProtected}
```

ImageLegend is an option for MakeSlicelImages, if set true a barlegend is added to the image. >>

```
Attributes[ImageLegend] = {Protected, ReadProtected}
```

ImageSize is an option that specifies the overall size of an image to display for an object. >>

```
Attributes[ImageSize] = {Protected}
```

MakeCheckPlot is an option for GetSlicePositions and if set true gives a plot of the slices locations. >>

```
Attributes[MakeCheckPlot] = {Protected, ReadProtected}
```

Method is an option for various algorithm-intensive functions that specifies what internal methods they should use. >>

```
Attributes[Method] = {Protected}
```

NormalizeIVIM is an option for IVIMplot. If True the signal at b=0 is 1. >>

```
Attributes[NormalizeIVIM] = {Protected, ReadProtected}
```

PeakNumber is an option of GetSlicePostitions and specifies how many slices per direction need to be found. >>

```
Attributes[PeakNumber] = {Protected, ReadProtected}
```

PlotColor is an option for GradientPlot can be any color or gradient color name. >>

```
Attributes[PlotColor] = {Protected, ReadProtected}
```

PlotRange is an option for graphics functions that specifies what range of coordinates to include in a plot. >>

```
Attributes[PlotRange] = {Protected, ReadProtected}
```

PlotSpace is an option for GradientPlot can be "bspace" or "qspace". >>

```
Attributes[PlotSpace] = {Protected, ReadProtected}
```

PositiveZ is an options for GradientPlot. If True all Gradients are displayed with a positive z direction. >>

```
Attributes[PositiveZ] = {Protected, ReadProtected}
```

SphereColor ListSpherePlor. Default value is Automatic, If a color is given this color will be used for all spheres. >>

```
Attributes[SphereColor] = {Protected, ReadProtected}
```

SphereSize is an option for GradientPlot and
ListSpherePlot. Sets the size of the spheres that represent the gradients. >>

```
Attributes[SphereSize] = {Protected, ReadProtected}
```

ProcessingTools

Functions

CorrectJoinSetMotion[[{dat1,dat2,...}, vox, over] motion corrects multiple sets with overlap. Over is the number of slices overlap between sets. A Translation registration is performed. >>

```
Attributes[CorrectJoinSetMotion] = {Protected, ReadProtected}
```

```
Options[CorrectJoinSetMotion] = {JoinSetSplit → True, PaddOverlap → 2}
```

```
SyntaxInformation[CorrectJoinSetMotion] = {ArgumentsPattern → {_, _, _}, OptionsPattern[]}
```

DataTransformation[data,vox,w] transforms a 3D dataset according to the affine transformation vector w

```
Attributes[DataTransformation] = {Protected, ReadProtected}
```

```
Options[DataTransformation] = {InterpolationOrder → 1}
```

```
SyntaxInformation[DataTransformation] = {ArgumentsPattern → {_, _, _}, OptionsPattern[]}
```

DatTot[{data1, data2, ..}, name, vox] calculates the parameter table containing the volume, mean, std and 95 CI for each of the diffusion parameters. >>

```
Attributes[DatTot] = {Protected, ReadProtected}
```

```
SyntaxInformation[DatTot] = {ArgumentsPattern → {_, _, _}}
```

DatTotXLS[{data1, data2, ..}, name, vox] is the same as
DatTot, but gives the parameters as strings for easy export to excel. >>

```
Attributes[DatTotXLS] = {Protected, ReadProtected}
```

```
SyntaxInformation[DatTotXLS] = {ArgumentsPattern → {_, _, _}}
```

ErrorPlot[data, xdata] plots a errorplot of the data where
the first dim of the data is the xrange which matches the xdata list.
ErrorPlot[data, xdata, range] similar with a given y range. >>

```
Attributes[ErrorPlot] = {Protected, ReadProtected}
```

```
Options[ErrorPlot] =  
{ColorValue → {GrayLevel[0], RGBColor[1, 0, 0]}, PlotLabel → , AxesLabel → , ImageSize → 300}
```

```
SyntaxInformation[ErrorPlot] = {ArgumentsPattern → {_, _, _}, OptionsPattern[]}
```

`FiberDensityMap[fiberPoin, dim, vox]` generates a fiber density map for the fiberPoin which are imported by LoadFiberTracts. The dimensions dim should be the dimensions of the tracked datasets van vox its voxel size. >>

```
Attributes[FiberDensityMap] = {Protected, ReadProtected}
```

```
Options[FiberDensityMap] = {SeedDensity → Automatic}
```

```
SyntaxInformation[FiberDensityMap] = {ArgumentsPattern → {_, _, _}, OptionsPattern[]}
```

`FiberLengths[fpoints,flines]` calculates the fiber length using the output from LoadFiberTracts.
`FiberLengths[{fpoints,flines}]` calculates the fiber length using the output from LoadFiberTracts. >>

```
Attributes[FiberLengths] = {Protected, ReadProtected}
```

```
SyntaxInformation[FiberLengths] = {ArgumentsPattern → {_, _}}
```

`FindOutliers[data]` finds the outliers of a list of data. >>

```
Attributes[FindOutliers] = {Protected, ReadProtected}
```

```
Options[FindOutliers] = {OutlierMethod → IQR, OutlierOutput → Mask,  
  OutlierIterations → 1, OutlierRange → 2, OutlierIncludeZero → True}
```

```
SyntaxInformation[FindOutliers] = {ArgumentsPattern → {_, _}, OptionsPattern[]}
```

`FitData[data,range]` converts the data into 100 bins within the +/- range around the mean. Function is used in ParameterFit. >>

```
Attributes[FitData] = {Protected, ReadProtected}
```

```
SyntaxInformation[FitData] = {ArgumentsPattern → {_, _}}
```

`GetMaskMeans[dat, mask, name]` calculates the mean, std, 5,50 and 95% CI form the given data for each of the given masks.
 Mask can be generated by SplitSegmentations. name is a string that is added to the header. >>

```
Attributes[GetMaskMeans] = {Protected, ReadProtected}
```

```
Options[GetMaskMeans] = {MeanMethod → SkewNormalDist}
```

```
SyntaxInformation[GetMaskMeans] = {ArgumentsPattern → {_, _, _}, OptionsPattern[]}
```

`Hist[data, range]` plots a probability density histogram of the data from xmin to xmax with a fitted (skew)normal distribution. Uses ParameterFit.

`Hist[data, range, label]` plots a probability density histogram of the data from xmin to xmax with a fitted (skew)normal distribution and label as x-axis label.

`Hist[{data1,..,data2,..}, {range1,range2,..}]` plots a probability density histogram of the data from xmin to xmax with a fitted (skew)normal distribution. Uses ParameterFit.

`Hist[{data1,data2,..}, {range1,range2,..}, {label1,label2,..}]` plots a probability density histogram of the data from xmin to xmax with a fitted (skew)normal distribution and label as x-axis label. >>

```
Attributes[Hist] = {Protected, ReadProtected}
```

```
Options[Hist] = {ColorValue →  
  {{GrayLevel[0], GrayLevel[1]}, RGBColor[1, 0, 0], RGBColor[0, 1, 0], RGBColor[0, 0, 1]},  
  Method → SkewNormal, PlotLabel → , AxesLabel → , ImageSize → 300}
```

```
SyntaxInformation[Hist] = {ArgumentsPattern → {_, _., OptionsPattern[]}}
```

Hist2[pars, range] plots a probability density histogram of the data over range with two fitted (skew)normal distribution. Uses ParameterFit2.
Hist2[pars, range, label] plots a probability density histogram of the data over range with two fitted (skew)normal distribution. Uses ParameterFit2. >>

```
Attributes[Hist2] = {Protected, ReadProtected}
```

```
Options[Hist2] = {Scaling → False}
```

```
SyntaxInformation[Hist2] = {ArgumentsPattern → {_, _, _, OptionsPattern[]}}
```

InvertDataset[data] inverts the data along the x y and z axes. In other words it is rotated around the origin such that $(x,y,z)=(-x,-y,-z)$ and $(0,0,0)=(0,0,0)$ >>

```
Attributes[InvertDataset] = {Protected, ReadProtected}
```

JoinSets[{dat1,dat2,...}, over] joins dat1, dat2, ... with over slices overlap.
JoinSets[{dat1,dat2,dat3...},{over1,over2,...}] joins dat1 and dat2 with over1 slices overlap, Joins dat2 and dat3 with over2 slices overlap and so on.
JoinSets[{dat1,dat2,...},{over,drop1,drop2,...}] joins dat1, dat2 with over slices overlap and drops drop1 slices for dat1 and drop2 from drop 2. >>

```
Attributes[JoinSets] = {Protected, ReadProtected}
```

```
Options[JoinSets] = {ReverseSets → True, ReverseData → True,  
  NormalizeSets → True, MotionCorrectSets → False, PaddOverlap → 2, JoinSetSplit → True}
```

```
SyntaxInformation[JoinSets] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

MeanRange[Range] calculates the medain (50%) and standard deviation (14% and 86%) range and reports it as a string. >>

```
Attributes[MeanRange] = {Protected, ReadProtected}
```

MeanStd[data] calculates the mean and standard deviation and reports it as a string. >>

```
Attributes[MeanStd] = {Protected, ReadProtected}
```

MedCouple[data] calculates the medcouple of a list of data. >>

```
Attributes[MedCouple] = {Protected, ReadProtected}
```


NumberTableForm[data] makes a right aligned table of the numbers with 3 decimal percision.
 NumberTableForm[data, n] makes a right aligned table of the numbers with n decimal percision. >>

Attributes[NumberTableForm] = {Protected, ReadProtected}

Options[NumberTableForm] = {TableMethod → NumberForm, TableAlignments → Automatic,
 TableDepth → ∞, TableDirections → Column, TableHeadings → None, TableSpacing → Automatic}

SyntaxInformation[NumberTableForm] = {ArgumentsPattern → {_, _., OptionsPattern[]}}

ParameterFit[data] fits a (skew)Normal probability density function to the data.

ParameterFit[{data1, data2,...}] fits a (skew)Normal
 probability density function to each of the datasets. Is used in Hist. >>

Attributes[ParameterFit] = {Protected, ReadProtected}

Options[ParameterFit] = {FitFunction → SkewNormal, FitOutput → Parameters, Method → Automatic}

SyntaxInformation[ParameterFit] = {ArgumentsPattern → {_, OptionsPattern[]}}

ParameterFit2[data] fits two skewNormal probaility density fucntions
 to the data. Assuming two compartments, one for fat and one for muscle. >>

Attributes[ParameterFit2] = {Protected, ReadProtected}

Options[ParameterFit2] = {FitOutput → BestFitParameters}

SyntaxInformation[ParameterFit2] = {ArgumentsPattern → {_, OptionsPattern[]}}

SetupDataStructure[dcmFolder] makes nii folders and generates nii
 files for a directory of dmc data where the data is structured per subject. >>

Attributes[SetupDataStructure] = {Protected, ReadProtected}

SmartMask[input] crates a smart mask of input, which is either
 the tensor or the tensor parameters calculated using ParameterCalc.

SmartMask[input, mask] crates a smart mask of input and used
 the mask as a prior selection of the input. >>

Attributes[SmartMask] = {Protected, ReadProtected}

Options[SmartMask] =
 {Strictness → 0.5, MaskCompartment → Muscle, SmartMethod → Continuous, SmartMaskOutput → mask}

SyntaxInformation[SmartMask] = {ArgumentsPattern → {_, _., OptionsPattern[]}}

SNRCalc[data,masksig,masknoise] calculates the Signal to noise ratio
 of the signal selected by masksig and the noise selected by masknoise. >>

Attributes[SNRCalc] = {Protected, ReadProtected}

SyntaxInformation[SNRCalc] = {ArgumentsPattern → {_, _, _}}

`SNRMapCalc[data1, noisemap]` calculates the signal to noise ratio of the data using $MN[data]/(1/\sqrt{\pi/2} \text{ sigma})$, where sigma is the local mean of the noise map assuming it is a rician distribution.

`SNRMapCalc[{data1, data2}]` calculates the signal to noise ratio from two identical images using $MN[data1, data2] / (.5 \text{ SQRT}[2] \text{ STDV}[data2 - data1])$.

`SNRMapCalc[{data1, .. dataN}]` calculates the signal to noise ratio of the data using MN/sigma where the mean signal MN is the average voxels value over all dynamics N and the sigma is the standard deviation over all dynamics N. >>

`Attributes[SNRMapCalc] = {Protected, ReadProtected}`

`Options[SNRMapCalc] = {OutputSNR → SNR, SmoothSNR → 2}`

`SyntaxInformation[SNRMapCalc] = {ArgumentsPattern → {_, __, __, OptionsPattern[]}}`

`SplitSets[data, Nsets, Nover]` splits the data in Nsets with Nover slices overlap. >>

`Attributes[SplitSets] = {Protected, ReadProtected}`

`Options[SplitSets] = {ReverseSets → False, ReverseData → True, PaddOverlap → 0}`

`SyntaxInformation[SplitSets] = {ArgumentsPattern → {_, __, __, OptionsPattern[]}}`

Options

`AxesLabel` is an option for graphics functions that specifies labels for axes. >>

`Attributes[AxesLabel] = {Protected}`

`ColorValue` is an option for `Hist` and `ErrorPlot`. Default {Black, Red}. >>

`Attributes[ColorValue] = {Protected, ReadProtected}`

`FitFunction` is an option for `ParameterFit`. Options are "Normal" or "SkewNormal". Indicates which function will be fitted. >>

`Attributes[FitFunction] = {Protected, ReadProtected}`

`FitOutput` is an option for `ParameterFit` and `ParameterFit2`. Option can be "Parameters", "Function" or "BestFitParameters". >>

`Attributes[FitOutput] = {Protected, ReadProtected}`

`ImageSize` is an option that specifies the overall size of an image to display for an object. >>

`Attributes[ImageSize] = {Protected}`

`InterpolationOrder` is an option for `Interpolation`, as well as `ListLinePlot`, `ListPlot3D`, `ListContourPlot`, and related functions, that specifies what order of interpolation to use. >>

```
Attributes[InterpolationOrder] = {Protected}
```

JoinSetSplit is an option ofr CorrectJoinSetMotion. If True
RegisterDataTransformSplit is used else RegisterDataTransform is used. >>

```
Attributes[JoinSetSplit] = {Protected, ReadProtected}
```

MaskCompartment is an option for SmartMask. Can be "Muscle" or "Fat". >>

```
Attributes[MaskCompartment] = {Protected, ReadProtected}
```

MeanMethod is an option for GetMaskMeans.
The option can be "NormalDist", "SkewNormalDist", or "Mean". >>

```
Attributes[MeanMethod] = {Protected, ReadProtected}
```

Method is an option for various algorithm-intensive
functions that specifies what internal methods they should use. >>

```
Attributes[Method] = {Protected}
```

MotionCorrectSets is an option for JoinSets. True motion
corrects the individual stacs before joining using CorrectJoinSetMotion. >>

```
Attributes[MotionCorrectSets] = {Protected, ReadProtected}
```

NormalizeSets is an option for JoinSets. True normalizes the individual stacs before joining. >>

```
Attributes[NormalizeSets] = {Protected, ReadProtected}
```

OutlierIncludeZero is an option for FindOutliers. If set
to True all values that are zero are ignored and considered outliers. >>

```
Attributes[OutlierIncludeZero] = {Protected, ReadProtected}
```

OutlierIterations is an option for FindOutliers.
Specifies how many iterations are used to find the outliers.
Each iteration the outliers are reevaluated on the data with the
previously found outliers already rejected. >>

```
Attributes[OutlierIterations] = {Protected, ReadProtected}
```

OutlierMethod is an option for FindOutliers. values can be "IQR", "SIQR" or "aIQR". "IRQ" is used for
normly distributed data, "SIQR" or "aIQR" are better for skewed distributions. >>

```
Attributes[OutlierMethod] = {Protected, ReadProtected}
```

OutlierOutput is an option for FindOutliers. If value is "Mask" it gives
a list of 1 for data and 0 for outliers. Else the output is {data, outliers}. >>

```
Attributes[OutlierOutput] = {Protected, ReadProtected}
```

OutlierRange is an option for FindOutliers. Specifies how many times the IQR is considered an outlier. >>

```
Attributes[OutlierRange] = {Protected, ReadProtected}
```

OutputSNR is an option for SNRMapCalc. >>

```
Attributes[OutputSNR] = {Protected, ReadProtected}
```

PaddOverlap is an option of CorrectJoinSetMotion
and JoinSets. it allows for extra motion in the z direction. >>

```
Attributes[PaddOverlap] = {Protected, ReadProtected}
```

PlotLabel is an option for graphics functions that specifies an overall label for a plot. >>

```
Attributes[PlotLabel] = {Protected}
```

ReverseData is an option for JoinSets. Reverses each
individual dataset given as input for the JoinSets function. True by default. >>

```
Attributes[ReverseData] = {Protected, ReadProtected}
```

ReverseSets is an option for JoinSets. Reverses the order of the datasets, False by default. >>

```
Attributes[ReverseSets] = {Protected, ReadProtected}
```

Scaling is an option for Hist2. Scales the individual fits of the fat and muscle compartment. >>

```
Attributes[Scaling] = {Protected, ReadProtected}
```

SeedDensity is an option for FiberDensityMap. The seedpoint spacing in mm. >>

```
Attributes[SeedDensity] = {Protected, ReadProtected}
```

SmartMaskOutput is an option for Smartmask. Can be set to "mask"
to output only the mask or "full" to also output the probability mask. >>

```
Attributes[SmartMaskOutput] = {Protected, ReadProtected}
```

SmartMethod is an option for SmartMask. This specifies
how the mask is generated. Can be "Continuous" or "Catagorical" >>

```
Attributes[SmartMethod] = {Protected, ReadProtected}
```

SmoothSNR is an option for SNRMapCalc. >>

```
Attributes[SmoothSNR] = {Protected, ReadProtected}
```

Strictness is an option for SmartMask value between 0 and 1. Higer values removes more data. >>

```
Attributes[Strictness] = {Protected, ReadProtected}
```

TableAlignments is an option for TableForm and MatrixForm which specifies how entries in each dimension should be aligned. >>

```
Attributes[TableAlignments] = {Protected}
```

TableDepth is an option for TableForm and MatrixForm which specifies the maximum number of levels to be printed in tabular or matrix format. >>

```
Attributes[TableDepth] = {Protected}
```

TableDirections is an option for TableForm and MatrixForm which specifies whether successive dimensions should be arranged as rows or columns. >>

```
Attributes[TableDirections] = {Protected}
```

TableHeadings is an option for TableForm and MatrixForm which gives the labels to be printed for entries in each dimension of a table or matrix. >>

```
Attributes[TableHeadings] = {Protected}
```

TableMethod is an option for NumberTableForm. It specifies which number form to uses. Values can be NumberForm, ScientificForm or EngineeringForm >>

```
Attributes[TableMethod] = {Protected, ReadProtected}
```

TableSpacing is an option for TableForm and MatrixForm which specifies how many spaces should be left between each successive row or column. >>

```
Attributes[TableSpacing] = {Protected}
```

RelaxometryTools

Functions

CalibrateEPGT2Fit[datan, times, angle]
calculates the Fat T2 relaxation that will be used in the EPGT2fit.

Outputs the fat T2 value. >>

```
Attributes[CalibrateEPGT2Fit] = {Protected, ReadProtected}
```

```
Options[CalibrateEPGT2Fit] =  
{EPGRelaxPars -> {{0, 100}, {20, 300}, {1400., 365.}}, EPGFitPoints -> 50}
```

```
SyntaxInformation[CalibrateEPGT2Fit] = {ArgumentsPattern -> {_, _, _}, OptionsPattern[]}
```

CreateT2Dictionary[{T1m, T1f}, {Necho, echoSpace}, angle] Creates a EPG signal dictionary used for EPGT2fit. Every dictionary that is defined is cached.

Output is {dictionary, vals} >>

Attributes[CreateT2Dictionary] = {Protected, ReadProtected}

Options[CreateT2Dictionary] = {DictB1Range → {0.5, 1.4, 0.01}, DictT2Range → {10., 70., 0.2}, DictT2fRange → {100., 200., 2.}}

SyntaxInformation[CreateT2Dictionary] = {ArgumentsPattern → {_, _, _}, OptionsPattern[]}

DictionaryMinSearch[dictionary, y] performs dictionary minimization of data y. dictionary is generated with CreateT2Dictionary.

Output is {{T2, B1}, fwfraction, residualError}. >>

Attributes[DictionaryMinSearch] = {Protected, ReadProtected}

SyntaxInformation[DictionaryMinSearch] = {ArgumentsPattern → {_, _, _}}

EPGSignal[{Necho, echoSpace}, {T1, T2}, {ex_angle, ref_angle}, B1] generates a EPG T2 curve with stimulated echos. T1, T2 and echoSpace are in ms, angel is in degree, B1 is between 0 and 1.

Output is the EPG Signal vector. >>

Attributes[EPGSignal] = {Protected, ReadProtected}

SyntaxInformation[EPGSignal] = {ArgumentsPattern → {_, _, _, _, _}}

EPGT2Fit[data, {Necho, detlaTE}, {exitation, refoucs}] fits the T2 based on Marty B et.al. Simultaneous muscle water T2 and fat fraction mapping using transverse relaxometry with stimulated echo compensation. Exitation and refocus are the RF pulse angles e.g. 90,180. They can also be a range of angeles over the slice profile as defined by GetSliceProfile.

Output is {{{T2map,B1Map},{wat, fat, fatMap}, residual},callibration} or {{T2map,B1Map},{wat, fat, fatMap}, residual} >>

Attributes[EPGT2Fit] = {Protected, ReadProtected}

Options[EPGT2Fit] = {DictB1Range → {0.5, 1.4, 0.01}, DictT2Range → {10., 60., 0.2}, DictT2fRange → {120., 170., 1.}, EPGRelaxPars → {1400., 365.}, EPGCalibrate → False, EPGFitPoints → 50, EPGMethod → dictionaryM, MonitorEPGFit → True, OutputCalibration → False, EPGSmoothB1 → True}

SyntaxInformation[EPGT2Fit] = {ArgumentsPattern → {_, _, _}, OptionsPattern[]}

NonLinearEPGFit[{vals, T2cons}, y] performs dictionary minimization of data y. vals = {{T1muscle, T1fat, T2fat}, {Necho, echoSpace, angle}}.

Output is {{T2, B1}, fwfraction, residualError}. >>

Attributes[NonLinearEPGFit] = {Protected, ReadProtected}

SyntaxInformation[NonLinearEPGFit] = {ArgumentsPattern → {_, _}}

T1Fit[data, TR] fits the T1 value to the data using a nonlinear method.

Output is {t1, apar, bpar} >>

Attributes[T1Fit] = {Protected, ReadProtected}

T1rhoFit[data, EchoTimes] fits the T1rho value to the data using linear or nonlinear methods.

Output is {S(0), T1rhomap}. >>

Attributes[T1rhoFit] = {Protected, ReadProtected}

Options[T1rhoFit] = {Method → Linear}

SyntaxInformation[T1rhoFit] = {ArgumentsPattern → {_, _, OptionsPattern[]}}

T2Fit[data, EchoTimes] fits the T2 value to the data using linear or nonlinear methods.

Output is {S(0), T2}. >>

Attributes[T2Fit] = {Protected, ReadProtected}

Options[T2Fit] = {Method → Linear}

SyntaxInformation[T2Fit] = {ArgumentsPattern → {_, _, OptionsPattern[]}}

TriExponentialT2Fit[data, EchoTimes] fits the T2 based on Azzabou N et.al. Validation of a generic approach to muscle water T2 determination at 3T in fat-infiltrated skeletal muscle. J. Magn. Reson. 2015. The fat T2 parameters are automatically estimated from the high signal voxels from the last echo.

Output is {{S(0), fatFraction, muscleFraction, T2map}, calibration} or {S(0), fatFraction, muscleFraction, T2map}. >>

Attributes[TriExponentialT2Fit] = {Protected, ReadProtected}

Options[TriExponentialT2Fit] = {OutputCalibration → False}

SyntaxInformation[TriExponentialT2Fit] = {ArgumentsPattern → {_, _, OptionsPattern[]}}

Options

DictB1Range is an option for CreateT2Dictionary and EPGT2Fit. It specifies the range and step of the B1 values in the dictionary {min, max, step}. >>

Attributes[DictB1Range] = {Protected, ReadProtected}

DictT2fRange is an option for CreateT2Dictionary and EPGT2Fit. It specifies the range and step of the T2 fat values in the dictionary {min, max, step} in ms. If a single value is given this fixed value is used as long as EPGCalibrate is False. >>

Attributes[DictT2fRange] = {Protected, ReadProtected}

DictT2Range is an option for CreateT2Dictionary and EPGT2Fit. It specifies the range and step of the T2 values in the dictionary {min, max, step} in ms. >>

Attributes[DictT2Range] = {Protected, ReadProtected}

EPGCalibrate is an option for EPGT2Fit. If set to True it does automatic calibration of the T2 fat relaxation time. >>

Attributes[EPGCalibrate] = {Protected, ReadProtected}

EPGFitPoints is an option for CalibrateEPGT2Fit and EPGT2Fit. Number of points is 200 by default. >>

Attributes[EPGFitPoints] = {Protected, ReadProtected}

EPGMethod is an option for EPGT2Fit. Values can be "NLLS", "dictionary" or "dictionaryM". >>

Attributes[EPGMethod] = {Protected, ReadProtected}

EPGRelaxPars is an option for EPGT2Fit. Needs to be {T1muscl, T1Fat, T2Fat} in ms, default is {1400,365,137}. >>

Attributes[EPGRelaxPars] = {Protected, ReadProtected}

EPGSmoothB1 is an option for EPGT2Fit. If set to True the B1 map of the fit will be smoothed after which the minimization is performed again but with a fixed B1. >>

Attributes[EPGSmoothB1] = {Protected, ReadProtected}

Method is an option for various algorithm-intensive functions that specifies what internal methods they should use. >>

Attributes[Method] = {Protected}

MonitorEPGFit show waitbar during EPGT2Fit. >>

Attributes[MonitorEPGFit] = {Protected, ReadProtected}

OutputCalibration is an option for EPGT2Fit and TriExponentialT2Fit. If true it outputs the calibration values. >>


```
Attributes[OutputCalibration] = {Protected, ReadProtected}
```

SimulationTools

Functions

AddNoise[data, noise] adds rician noise to the data with a given sigma or SNR value. >>

```
Attributes[AddNoise] = {Protected, ReadProtected}
```

```
Options[AddNoise] = {NoiseSize → Sigma}
```

```
SyntaxInformation[AddNoise] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

BlochSeries[vectorIn, deltata, freqRange, B1] performs a Bloch simulation of an RF pulse. >>

```
Attributes[BlochSeries] = {Protected, ReadProtected}
```

CalculateGfactor[factors, sensitivity, Wmat] calculates a gfactor for given sensitivity maps and noise correlation W. given the sense factors which is a list of three integers. >>

```
Attributes[CalculateGfactor] = {Protected, ReadProtected}
```

```
Options[CalculateGfactor] = {GRegularization → 0.}
```

```
SyntaxInformation[CalculateGfactor] = {ArgumentsPattern → {_, _, _, _, OptionsPattern[]}}
```

CreateDiffData[sig, eig, bvec, gradients, dim] creates a DTI datasets of dimensions dim with sig as unweighted signal S0 and bvec and gradients. eig can be {l1, l2, l3}, {{l1, l2, l3}, {e1, e2, e3}}, {{l1, l2, l3}, "Random"}, {{l1, l2, l3}, "RandomZ"} or {{l1, l2, l3}, "OrtRandom"}.

Uses Tensor internally. >>

```
Attributes[CreateDiffData] = {Protected, ReadProtected}
```

```
SyntaxInformation[CreateDiffData] = {ArgumentsPattern → {_, _, _, _, _}}
```

GetPulseProfile[excitation, refocus] gives

the pulse angle profiles for the excitation and refocussing pulses.

a pulse is defined as {"name", flipangle, {G_strnth, Dur, BW}}.

GetPulseProfile[{"name", flipangle,

{G_strnth, Dur, BW}}] gives detailed slice profile information of one pulse.

output is {ex_angle_profile, ref_angle_profile, {plots}}.

output for single pulse is {{distance, Mt, Mz, Mx, My, ang, phase}, plots} >>

```
Attributes[GetPulseProfile] = {Protected, ReadProtected}
```

```
Options[GetPulseProfile] =
```

```
{MagnetizationVector → {0, 0, 1}, SliceRange → 12, SliceRangeSamples → 25}
```

```
SyntaxInformation[GetPulseProfile] = {ArgumentsPattern → {_, _, OptionsPattern[]}}
```

`GfactorSimulation[sensitivity, Wmat, {dir,sense}]` calculates the gfactormaps for given sensitivity maps and noise corraltion W in one direction. The sensefactors are a list of integers in a given direction: "LR", "FH", or "AP".
`GfactorSimulation[sensitivity, Wmat, {dir1,sense1}, {dir2,sense2}]` calculates the gfactormaps for given sensitivity maps and noise corraltion W in two directions. >>

`Attributes[GfactorSimulation] = {Protected, ReadProtected}`

`Options[GfactorSimulation] = {GRegularization → 0., GOutput → Grid}`

`SyntaxInformation[GfactorSimulation] = {ArgumentsPattern → {_, _, _, _., _., OptionsPattern[]}}`

`PlotSimulation[pars, xval, true, label, color]` plots the pars (output form Parameters). Using label as PlotLabel and xval as x axis Thics.tr are the true parameter values. color are the color used for the plot. >>

`Attributes[PlotSimulation] = {Protected, ReadProtected}`

`Options[PlotSimulation] = {PlotRange → {{0, 3}, {0, 3}, {0, 3}, {0, 3}, {0, 1}}}`

`SyntaxInformation[PlotSimulation] = {ArgumentsPattern → {_, _, _, _, _, OptionsPattern[]}}`

`PlotSimulationAngle[par, xdata, label, col]` plots pars (output from Anlge Parameters). >>

`Attributes[PlotSimulationAngle] = {Protected, ReadProtected}`

`Options[PlotSimulationAngle] = {PlotRange → {0, 90}}`

`PlotSimulationAngleHist[pars, label, xdata]` plots pars (output from Anlge Parameters). >>

`Attributes[PlotSimulationAngleHist] = {Protected, ReadProtected}`

`SyntaxInformation[PlotSimulationAngleHist] = {ArgumentsPattern → {_, _, _, _, _., OptionsPattern[]}}`

`PlotSimulationHist[pars, label, xdata, tr]` plots the pars (output form Parameters). Using label as plotlabel and xdata as x axis label. tr are the true parameter values. >>

`Attributes[PlotSimulationHist] = {Protected, ReadProtected}`

`SyntaxInformation[PlotSimulationHist] = {ArgumentsPattern → {_, _, _, _}}`

`PlotSimulationVec[tens, xdata, label]` plots the eigenvectors from simulated tensors. >>

`Attributes[PlotSimulationVec] = {Protected, ReadProtected}`

`Options[PlotSimulationVec] = {SortVecs → True}`

`SyntaxInformation[PlotSimulationVec] = {ArgumentsPattern → {_, _, _, OptionsPattern[]}}`

`Pulses[name]` gives the pulse shape of some predefinec Philips pulse shapes. >>

`Attributes[Pulses] = {Protected, ReadProtected}`

Signal[par,TR,TE] calculates the MRI signal at a given TR and TE. Par is defined as {pd, T1, T2}. >>

Attributes[Signal] = {Protected, ReadProtected}

SyntaxInformation[Signal] = {ArgumentsPattern → {_, _, _}}

SimAngleParameters[tens,vec] calculates the diffusion eigenvectors for tens compared to the true values vec. The output can be used in PlotSimulationAngleHist and PlotSimulationAngle. >>

Attributes[SimAngleParameters] = {Protected, ReadProtected}

SyntaxInformation[SimAngleParameters] = {ArgumentsPattern → {_, _}}

SimParameters[tens] calculates the diffusion parameters for tens. The output can be used in PlotSimulationHist and PlotSimulation. >>

Attributes[SimParameters] = {Protected, ReadProtected}

Options[SimParameters] = {Reject → False}

SyntaxInformation[SimParameters] = {ArgumentsPattern → {_, OptionsPattern[]}}

SimulateSliceEPG[exitation, refocus, {{T1, T2}, {Necho, echoSp}, b1}] gives a simulated slice profile and EPG signal plot. exitation and refocus are generated by GetPulseProfile. >>

Attributes[SimulateSliceEPG] = {Protected, ReadProtected}

Options[SimulateSliceEPG] = {ReportFits → False}

SyntaxInformation[SimulateSliceEPG] = {ArgumentsPattern → {_, _, _, OptionsPattern[]}}

Tensor[{l1, l2, l3}] creates a diffusion tensor with vectors {{0,0,1},{0,1,0},{1,0,0}} and eigenvalues {l1, l2, l3}.
 Tensor[{l1, l2, l3}, {e1, e2, e3}] creates a diffusion tensor with vectors {e1, e2, e3} and eigenvalues {l1, l2, l3}.
 Tensor[{l1, l2, l3}, "Random"] creates a diffusion tensor with random orthogonal eigenvectors {e1, e2, e3} and eigenvalues {l1, l2, l3}.
 Tensor[{l1, l2, l3}, "RandomZ"] creates a diffusion tensor with random orthogonal eigenvectors {{1,0,0}, e2, e3} with random eigenvectors and eigenvalues {l1, l2, l3}.
 Tensor[{l1, l2, l3}, "OrtRandom"] creates a diffusion tensor with random orthogonal eigenvectors {{1,0,0},{0,1,0},{0,0,1}} and eigenvalues {l1, l2, l3}. >>

Attributes[Tensor] = {Protected, ReadProtected}

Options[Tensor] = {TensOutput → Vector}

SyntaxInformation[Tensor] = {ArgumentsPattern → {_, .., OptionsPattern[]}}

Options

GOutput is an option for GfactorSimulation. can be "Grid" or "List". >>

```
Attributes[GOutput] = {Protected, ReadProtected}
```

GRegularization is an option for CalculateGfactor and GfactorSimulation. >>

```
Attributes[GRegularization] = {Protected, ReadProtected}
```

MagnetizationVector is an option for GetPulseProfile.

It defines the start magnetization vector for the bloch simulation. >>

```
Attributes[MagnetizationVector] = {Protected, ReadProtected}
```

NoiseSize is an option for AddNoise. Values can be

"Sigma", then the noise sigma is given or "SNR", then the SNR is given. >>

```
Attributes[NoiseSize] = {Protected, ReadProtected}
```

PlotRange is an option for graphics functions

that specifies what range of coordinates to include in a plot. >>

```
Attributes[PlotRange] = {Protected, ReadProtected}
```

Reject is an option for EigenvalCalc. If True

then voxels with negative eigenvalues are rejected and set to 0. >>

```
Attributes[Reject] = {Protected, ReadProtected}
```

ReportFits is an option for SimulateSliceEPG. If True it also reports the fit values >>

```
Attributes[ReportFits] = {Protected, ReadProtected}
```

SliceRange is an option for GetPulseProfile. It specifies over which range

the slice profile is generated (in mm). the total profile is 2xSliceRange. >>

```
Attributes[SliceRange] = {Protected, ReadProtected}
```

SliceRangeSamples is an option for GetPulseProfile.

defines how many samples are used to generate half a puls profile. >>

```
Attributes[SliceRangeSamples] = {Protected, ReadProtected}
```

SortVecs is an option for PlotSimulationVec. >>

```
Attributes[SortVecs] = {Protected, ReadProtected}
```

TensOutput is an option for Tensor. Values can be "Vector" or "Matrix". >>

```
Attributes[TensOutput] = {Protected, ReadProtected}
```

TensorTools

Functions

ADCCalc[eigenvalues] calculates the ADC from the given eigenvalues. >>

```
Attributes[ADCCalc] = {Protected, ReadProtected}
```

```
SyntaxInformation[ADCCalc] = {ArgumentsPattern -> { _ } }
```

AngleCalc[data, vector] calculates the angel between the vector and the data. Data should be an array of dimensions {xxx,3}. >>

```
Attributes[AngleCalc] = {Protected, ReadProtected}
```

```
Options[AngleCalc] = {Distribution -> 0-180}
```

```
SyntaxInformation[AngleCalc] = {ArgumentsPattern -> { _, _, OptionsPattern[] } }
```

AngleMap[data] calculates the zennith and azimuth angles of a 3D dataset (z,x,y,3) containing vectors relative to the slice direction. >>

```
Attributes[AngleMap] = {Protected, ReadProtected}
```

```
SyntaxInformation[AngleMap] = {ArgumentsPattern -> { _ } }
```

ColorFAPlot[tenor] create a color coded FA map from the tensor for l1, l2 and l3. >>

```
Attributes[ColorFAPlot] = {Protected, ReadProtected}
```

```
SyntaxInformation[ColorFAPlot] = {ArgumentsPattern -> { _ } }
```

ConcatenateDiffusionData[{data1, ..., dataN}, {grad1, ..., gradN}, {bval, ..., bvalN}, {vox, ..., voxN}] concatenates the diffusion data sets.
ConcatenateDiffusionData[{data1, ..., dataN}, {grad1, ..., gradN}, {bval, ..., bvalN}, {vox, ..., voxN}] concatenates the diffusion data sets. >>

```
Attributes[ConcatenateDiffusionData] = {Protected, ReadProtected}
```

```
SyntaxInformation[ConcatenateDiffusionData] = {ArgumentsPattern -> { _, _, _, _ } }
```

Correct[data, phase, shiftpar] corrects the dataset data using the phasemap and the shiftpar and interpolation order 1.
Correct[data, phase, shiftpar, int] corrects the dataset data using the phasemap and the shiftpar and interpolation order int. >>

```
Attributes[Correct] = {Protected, ReadProtected}
```

```
SyntaxInformation[Correct] = {ArgumentsPattern -> { _, _, _, _ } }
```

Deriv[disp, vox] calculates the derivative of the displacement along the three main axes. disp is the displacement field, vox is the voxel size.
Deriv[disp, vox, mask] calculates the derivative of the displacement along the three main axes. Sharp edges between the background en disp are solved by the mask. mask is a mask delining the edge of the displacement field. >>

```
Attributes[Deriv] = {Protected, ReadProtected}
```

```
SyntaxInformation[Deriv] = {ArgumentsPattern → {_, _, _}}
```

DriftCorrect[data, bval] dirft corrects the data using
the signals of the lowest bvalue that has 6 or more unique volumes.
For the function to work optimal it is best to have these volumes evenly spread
througout thet data and for the first and last volume to have this low bvalue. >>

```
Attributes[DriftCorrect] = {Protected, ReadProtected}
```

```
Options[DriftCorrect] = {NormalizeSignal → True, UseMask → True}
```

```
SyntaxInformation[DriftCorrect] = {ArgumentsPattern → {_, _, _}, OptionsPattern[]}
```

ECalc[eigenvalues] caculates the e from the given eigenvalues. >>

```
Attributes[ECalc] = {Protected, ReadProtected}
```

```
Options[ECalc] = {MonitorCalc → True}
```

```
SyntaxInformation[ECalc] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

EigensysCalc[tensor] caculates the eigensystem for the given tensor. >>

```
Attributes[EigensysCalc] = {Protected, ReadProtected}
```

```
Options[EigensysCalc] = {MonitorCalc → False}
```

```
SyntaxInformation[EigensysCalc] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

EigenvalCalc[tensor] caculates the eigenvalues for the given tensor. >>

```
Attributes[EigenvalCalc] = {Protected, ReadProtected}
```

```
Options[EigenvalCalc] = {MonitorCalc → True, RejectMap → False, Reject → True}
```

```
SyntaxInformation[EigenvalCalc] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

EigenvecCalc[tensor] caculates the eigenvectors for the given tensor. >>

```
Attributes[EigenvecCalc] = {Protected, ReadProtected}
```

```
Options[EigenvecCalc] = {MonitorCalc → False}
```

```
SyntaxInformation[EigenvecCalc] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

FACalc[eigenvalues] caculates the FA from the given eigenvalues. >>

```
Attributes[FACalc] = {Protected, ReadProtected}
```

```
SyntaxInformation[FACalc] = {ArgumentsPattern → {_}}
```

ParameterCalc[tensor] caculates the eigenvalues and MD and FA from the given tensor. >>

```
Attributes[ParameterCalc] = {Protected, ReadProtected}
```

```
Options[ParameterCalc] = {Reject → True, MonitorCalc → False}
```

```
SyntaxInformation[ParameterCalc] = {ArgumentsPattern → {_, OptionsPattern[]}}
```

RemoveIsoImages[data, grad, bval] Removes the ISO images from the philips scanner from the data. ISO images have $g=\{0,0,0\}$ and $b>0$. >>

```
Attributes[RemoveIsoImages] = {Protected, ReadProtected}
```

```
SyntaxInformation[RemoveIsoImages] = {ArgumentsPattern → {_, _, _}}
```

ResidualCalc[DTI,{tensor,S0},gradients,bvector] calculates the tensor residuals for the given dataset.

ResidualCalc[DTI,{tensor,S0},outlier,gradients,bvector] calculates the tensor residuals for the given dataset taking in account the outliers.

ResidualCalc[DTI,{tensor,S0},bmat] calculates the tensor residuals for the given dataset.

ResidualCalc[DTI,{tensor,S0},outlier,bmat] calculates the tensor residuals for the given dataset taking in account the outliers.

ResidualCalc[DTI,tensor,gradients,bvector] calculates the tensor residuals for the given dataset. Tensor must contain Log[S0].

ResidualCalc[DTI,tensor,outlier,gradients,bvector] calculates the tensor residuals for the given dataset taking in account the outliers. Tensor must contain Log[S0].

ResidualCalc[DTI,tensor,bmat] calculates the tensor residuals for the given dataset. Tensor must contain Log[S0].

ResidualCalc[DTI,tensor,outlier,bmat] calculates the tensor residuals for the given dataset taking in account the outliers. Tensor must contain Log[S0]. >>

```
Attributes[ResidualCalc] = {Protected, ReadProtected}
```

```
Options[ResidualCalc] = {MeanRes → All}
```

```
SyntaxInformation[ResidualCalc] = {ArgumentsPattern → {_, _, _, _, OptionsPattern[]}}
```

SigmaCalc[DTI,grad,bvec] calculates the noise sigma based on the tensor residual, using a blur factor of 10.

SigmaCalc[DTI,tens,grad,bvec] calculates the noise sigma based on the tensor residual, using a blur factor of 10.

SigmaCalc[DTI,grad,bvec,blur] calculates the noise sigma based on the tensor residual, If blur is 1 ther is no blurring.

SigmaCalc[DTI,tens,grad,bvec,blur] calculates the noise sigma based on the tensor residual. If blur is 1 ther is no blurring. >>

```
Attributes[SigmaCalc] = {Protected, ReadProtected}
```

```
Options[SigmaCalc] = {FilterShape → Median}
```

```
SyntaxInformation[SigmaCalc] = {ArgumentsPattern → {_, _, _, _, OptionsPattern[]}}
```

SortDiffusionData[data, grad, bval] sorts the diffusion datasets grad and bval for magnitude of bvalue. >>

Attributes[SortDiffusionData] = {Protected, ReadProtected}

SyntaxInformation[SortDiffusionData] = {ArgumentsPattern → {_, _, _}}

TensorCalc[data, gradients, bvalue] calculates the diffusion tensor for the given dataset. Allows for one unweighted image and one b value. Gradient directions must be in the form {{x1,y1,z1}, ..., {xn,yn,zn}} without the unweighted gradient direction. bvalue is a single number indicating the b-value used. **TensorCalc**[data, gradients, bvec] calculates the diffusion tensor for the given dataset. allows for multiple unweighted images and multiple bvalues. allows for different tensor fitting methods. gradient directions must be in the form {{x1,y1,z1}, ..., {xn,yn,zn}} with the unweighted direction as {0,0,0}. bvec the bvector, with a bvalue defined for each gradient direction. b value for unweighted images is 0. **TensorCalc**[data, bmatix] calculates the diffusion tensor for the given dataset. allows for multiple unweighted images and multiple bvalues. bmat is the bmatrix which can be generated using Bmatrix. >>

Attributes[TensorCalc] = {Protected, ReadProtected}

Options[TensorCalc] = {MonitorCalc → True, Method → iWLLS, FullOutput → True, RobustFit → True, Parallelize → True, RobustFitParameters → {0.0001, 6}}

SyntaxInformation[TensorCalc] = {ArgumentsPattern → {_, _, _, OptionsPattern[]}}

TensorCorrect[tensor, phase, shift, vox] corrects the tensor based on B0 field map. Can perform both translation and rotation of tensor. >>

Attributes[TensorCorrect] = {Protected, ReadProtected}

Options[TensorCorrect] = {RotationCorrect → False}

SyntaxInformation[TensorCorrect] = {ArgumentsPattern → {_, _, _, _, OptionsPattern[]}}

Options

Distribution is an option for AngleCalc. values can be "0-180", "0-90" and "-90-90". >>

Attributes[Distribution] = {Protected, ReadProtected}

FilterShape is an option for SigmaCalc. Can be "Gaussian" or "Median". >>

Attributes[FilterShape] = {Protected, ReadProtected}

FullOutput is an option for TensorCalc when using bvector. When True also the S0 is given as output. >>

Attributes[FullOutput] = {Protected, ReadProtected}

MeanRes is an option for ResidualCalc.

When True the root mean square of the residual is calculated. >>

```
Attributes[MeanRes] = {Protected, ReadProtected}
```

Method is an option for various algorithm-intensive

functions that specifies what internal methods they should use. >>

```
Attributes[Method] = {Protected}
```

MonitorCalc is an option for all Calc functions. When true the processes of the calculation is shown. >>

```
Attributes[MonitorCalc] = {Protected, ReadProtected}
```

NormalizeSignal is an option for DriftCorrect. >>

```
Attributes[NormalizeSignal] = {Protected, ReadProtected}
```

Parallelize[*expr*] evaluates *expr* using automatic parallelization. >>

```
Attributes[Parallelize] = {HoldFirst, Protected}
```

```
Parallelize[Parallel`Kernels`Private`args$___] :=  
(Parallel`Protected`doAutoLaunch[TrueQ[Parallel`Static`$enableLaunchFeedback]];  
Parallelize[Parallel`Kernels`Private`args$])
```

```
Options[Parallelize] = {Method -> Automatic, DistributedContexts -> $Context}
```

Reject is an option for EigenvalCalc. If True

then voxels with negative eigenvalues are rejected and set to 0. >>

```
Attributes[Reject] = {Protected, ReadProtected}
```

RejectMap is an option for EigenvalCalc. If Reject is True and RejectMap is True

both the eigenvalues as well as a map showing rejected values is returned. >>

```
Attributes[RejectMap] = {Protected, ReadProtected}
```

RobustFit is an option for TensorCalc. If true outliers will be rejected in the fit, only works with WLLS.

If FullOutput is given the outlier map is given. >>

```
Attributes[RobustFit] = {Protected, ReadProtected}
```

RobustFitParameters is an option for TensorCalc. gives the threshold for

stopping the iterations and the kappa for the outlier margin, {tr,kappa}. >>

```
Attributes[RobustFitParameters] = {Protected, ReadProtected}
```

RotationCorrect is an option for TensorCorrect. Default is False. Is a tensor is deformed

setting to True also the shear is accounted for by local rotation of the tensor >>

```
Attributes[RotationCorrect] = {Protected, ReadProtected}
```

UseMask is a function for MeanSignal and DriftCorrect >>

```
Attributes[UseMask] = {Protected, ReadProtected}
```

VisteTools

Functions

DatRead[file] imports data from file (dtitool *.dat format) as binary data using Real32 format. >>

```
Attributes[DatRead] = {Protected, ReadProtected}
```

```
SyntaxInformation[DatRead] = {ArgumentsPattern -> { _ } }
```

DatWrite[file, data] exports data to *.dat file as binary data using Real32 format. >>

```
Attributes[DatWrite] = {Protected, ReadProtected}
```

```
SyntaxInformation[DatWrite] = {ArgumentsPattern -> { _, _ } }
```

DTItoolExp[tensor, voxsize] exports tensor to {XX.dat, YY.dat, ZZ.dat, XY.dat, XZ.dat, YZ.dat} and uses XX.dat as background and generates corresponding *.dti files.

DTItoolExp[tensor, voxsize, folder] exports tensor to {XX.dat, YY.dat, ZZ.dat, XY.dat, XZ.dat, YZ.dat} to the given folder and uses XX.dat as background and generates corresponding *.dti files.

DTItoolExp[tensor, voxsize, folder, add] exports tensor to {XX.dat, YY.dat, ZZ.dat, XY.dat, XZ.dat, YZ.dat} to the given folder and uses XX.dat as background and generates corresponding *.dti files adds – add to the filenames.

DTItoolExp[back, tensor, voxsize] exports background to back.dat and tensor to {XX.dat, YY.dat, ZZ.dat, XY.dat, XZ.dat, YZ.dat} and generates corresponding *.dti files.

DTItoolExp[back, tensor, voxsize, folder] exports background to back.dat and tensor to {XX.dat, YY.dat, ZZ.dat, XY.dat, XZ.dat, YZ.dat} to the given folder and generates corresponding *.dti files.

DTItoolExp[back, tensor, voxsize, folder, add] exports background to back.dat and tensor to {XX.dat, YY.dat, ZZ.dat, XY.dat, XZ.dat, YZ.dat} to the given folder and generates corresponding *.dti files and adds – add to the filenames. >>

```
Attributes[DTItoolExp] = {Protected, ReadProtected}
```

```
SyntaxInformation[DTItoolExp] = {ArgumentsPattern -> { _, _, _, _, _ } }
```

DTItoolExpFile[file, background, add, voxsize] exports a *.dti text file. >>

```
Attributes[DTItoolExpFile] = {Protected, ReadProtected}
```

```
SyntaxInformation[DTItoolExpFile] = {ArgumentsPattern -> { _, _, _, _ } }
```

DTItoolExpInd[data, file] exports a 3D array

data to the file filename DTItool format (*.dat) using DatWrite.

DTItoolExpInd[data, file, folder] exports data to given file and folder.

DTItoolExpInd[data, file, folder, add]

exports data to given file and folder and adds –add to the filename. >>

